

Product datasheet for **RC231346**

CNOT4 (NM_001190847) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CNOT4 (NM_001190847) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CNOT4
Synonyms:	CLONE243; NOT4; NOT4H
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RC231346 representing NM_001190847
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCTCGCAGTCTGATGCGAAGGAAGACCCTGTGGAGTGCCCTCTTTCATGGAGCCCTTGGAGATAG
ATGATATCAACTTTTCCCTTGACCTGTGGCTACCAGATTTGCCGATTTGTTGGCATCGAATTCGCAC
TGATGAAAATGGGCTTTGTCTGATGTAGAAAGCCATATCCAGAAGACCCAGCAGTTTATAAAACCACTC
TCCAGGAAGAGCTGCAAAGGATAAAGAATGAGAAAAACAGAAACAAAATGAGAGAAAACAGAAAATAT
CAGAAAATCGCAAACATTTGGCTAGTGTACGTGTCGTACAAAAAACCTCGTCTTTGTTGTAGGTTTATC
TCAGCGCTAGCAGACCCAGAGGTTTTAAAAACGACCAGAATATTTTGGGAAGTTTGGTAAAAACATAAA
GTTGTCATCAATAATAGCACATCATATGCAGGCTCACAGGGTCCAAGTCCAGTGCTTATGTAACCTATA
TCCGGTCAGAAGACGCTCTCAGAGCCATACAGTGTGTCAACAATGTGGTAGTAGTGGCAGAACACTTAA
GGCATCTCTAGGTACAACAAAATACTGCAGTTACTTCTTAAAGAATATGCAGTGTCCAAAACCTGACTGC
ATGATCTTCATGAATTTGGGGATGAGGCGGCCAGCTTACAAAAGAGGAAATGCAGGCGGGTAAACACC
AAGAATATGAACAGAAGCTACTTCAAGAATTATATAAAATTAATCCCAATTTTCTTCAGCTATCTACGGG
TTCAGTTGATAAAAAAAGAACAAGTGACACCACTGCAGAGGTACGATACCCCCATTGACAAAACCTTCA
GATTCTCTCAGTATAGGGAACGGTGATAATCCCAGCAGATATCTAACAGTGTACGCCTTACCACCAC
CTGGTTTGTCAAATCCAATCCAGTCAATCCCATCAGTTCATCCAATCACAGTGCACGGTCCCCTTTTGA
AGGGGCAGTAACAGAGTCACAGTCGTTATTCTCAGACAATTTTCGCCATCCCAACCCTATCCCAAGTGGG
CTTCTCTCTTTCCCAGCTCCCCACAGACATCCAGTACTGGCCTACAGCACCAGAACCACAGAGCCTCT
TCACATCAGAAACAATCCAGTATCATCCTCTACAGACTGGCAAGCAGCTTTTGGCTTTGGTTCTTCTAA
ACAACCAGAGGATGACTTGGGTTTTGATCCCTTCGATGTCACTCGAAAAGCCTTAGCAGACCTGATTGAG
AAGGAAGTGTCCGTTCAAGACCAACCTTCCCTTTCGCCACATCTTTCAGAAGTCCCTTTCACACACTA
CAACCGCCAAAGGTCCAGGCTCTGGATTCTGCATCCTGCTGCAGCTACAAATGCCAATTCTCTCAATAG
TACCTTTTCAGTCTTGGCCAGAGGTTCCCTCAATTTTCAGCAGCACCGAGCGGTTATAATTCAATCAGT
TTTCCAGGCCAGGACCCGCTATCCTTGGATGGCCTTTCACGCAATAGCATCATGCATTGAACCACA
CAGCAAACCCACCTCAAATAGTAATTTCTTGGACTTGAATCTCCCGCCACAGCACAACACAGGTCTGGG
AGGGATCCCTGTAGCAGGTATTCAGCGTCTTCAGGAAACAGTTTAGACTCTCTCAAGATGACAATCCT
CCACTGGCTAAAATCCCTTCCAGCCCTCACAGAGATGGACGGCCCCAGCGCTGCTCCATCACAGACCC
ACCACAGCGCCCCCTTACGACACAGATCCCGCTGCACAGAGCCAGTTGGAATCCCTACCCTCCTCCTTC
AAACCTTCCAGTCTTCCACTCCCCACCCCAAGGCTTTCAGACAGCCTTTCAGACCCCCAGAAAACCCCC
ACAGATTTACTACAGAGTTCAACTGGACCGCCAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231346 representing NM_001190847
Red=Cloning site Green=Tags(s)


```
MSRSPDAKEDPVECPLCMEPLEIDDINFFPCTCGYQICRFCHWRIRTDENGLCPACRKPYPEDPAVYKPL
SQEELQRIKNEKKQKQNERKQKISENRKHLASVRVVQKNLVFVVLGSQRLADPEVLKRPEYFGKFGKIHK
VVINNSTSYAGSQGPSASAYVTYIRSEDALRAIQCVNNVVVDGRTLKASLGTTKYCSYFLKNMQCPKPCD
MYLHELGDDEAASFKEEMQAGKHQEQYEQKLLQELYKLNPNFLQLSTGSDVKNKNKVTPLQRYDTPIDKPS
DSLISGNGDNSQQISNSDTPSPPPGLSKSNPVIPISSSNHSARSPFEGAVTESQSLSFDNFRHNPPIPSG
LPPFPSSPQTSSDWPTAPEPQSFLTSETIPVSSSTDWQAAFQFGSSKQPEDDLGFPDVTTRKALADLIE
KELSVQDQPSLSPTLQNSSHTTTAKGPGSGFLHPAAATNANSLNSTFVLVQRFPQFQQHRAVYNSFS
FPGQAARYPWMAFPRNSIMHLNHTANPTSNSNFLDLNLPQHNTLGGIPVAGIPASSGNSLDSLQDDNP
PHWLKSLQALTEMDGPSAAPSQTHHSAPFSTQIPLHRASWNPYPSPSNPSSFHSPPPGFQTAFRPPSKTP
TDLLQSSTLDRH
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8071_b04.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001190847

ORF Size: 1926 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001190847.2](#)

RefSeq ORF: 1929 bp

Locus ID: 4850

UniProt ID: [O95628](#)

Cytogenetics: 7q33


Protein Families: Druggable Genome, Transcription Factors

Protein Pathways: RNA degradation

MW: 71.2 kDa

Gene Summary: The protein encoded by this gene is a subunit of the CCR4-NOT complex, a global transcriptional regulator. The encoded protein interacts with CNOT1 and has E3 ubiquitin ligase activity. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jul 2010]

Product images:



Circular map for RC231346