

Product datasheet for **RC231363L3V**

FAM59B (GAREM2) (NM_001191033) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FAM59B (GAREM2) (NM_001191033) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FAM59B
Synonyms:	FAM59B; GAREML
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001191033
ORF Size:	1992 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231363).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001191033.1 , NP_001177962.1
RefSeq ORF:	1995 bp
Locus ID:	150946
UniProt ID:	Q75VX8
Cytogenetics:	2p23.3
MW:	71.8 kDa
Gene Summary:	Probable adapter protein that may provide a link between cell surface epidermal growth factor receptor and the MAPK/ERK signaling pathway.[UniProtKB/Swiss-Prot Function]



[View online »](#)