

Product datasheet for **RC231390**

TAF6 (NM_001190415) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TAF6 (NM_001190415) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TAF6
Synonyms:	ALYUS; MGC:8964; TAF(II)70; TAF(II)80; TAF2E; TAFII-70; TAFII-80; TAFII70; TAFII80; TAFII85
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC231390 representing NM_001190415
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGACAAAGCGCCGTCTACAAACAATCACCGCTTGCCTTCAACTTCGGGAAGGGTCTCCGTCTCTCCACC
 GGGGGTTCATCCTTCCAGGGAGGAGAAGAGGGACTCCAGAATGGCTGAGGAGAAGAAGCTGAAGCTTAG
 CAACACTGTGCTGCCCTCGGAGTCCATGAAGGTGGTGGCTGAATCCATGGGCATCGCCAGATTCAGGAG
 GAGACCTGCCAGCTGCTAACGGATGAGGTGAGCTACCGCATCAAAGAGATCGCACAGGATGCCTTGAAGT
 TCATGCACATGGGGAAGCGGCAGAAGCTCACCACAGTGACATTGACTACGCCTTGAAGCTAAAGAATGT
 CGAGCCACTCTATGGCTTCCACGCCAGGATTCATTCTTTCCGCTTCGCCTCTGGTGGGGCCGGGAG
 CTTTACTTCTATGAGGAGAAGGAGTTGATCTGAGCGACATCATCAATACCCCTCTGCCCGGGTGGCCC
 TGGACGTCTGCCTCAAAGCTCATTGGCTGAGCATCGAGGGCTGCCAGCCAGCTATCCCGAGAACCAGCC
 CCCAGCTCCCAAAGAGCAACAGAAGGCTGAAGCCACAGAACCCTGAAGTCAGCCAAGCCAGGCCAGGAG
 GAAGACGGACCCCTGAAGGGCAAAGGTCAAGGGGCCACCACAGCCGACGGCAAAGGGAAAGAGAAGAAGG
 CGCCGCCCTTGCTGGAGGGGGCCCCCTTGCAGCTGAAGCCCCGGAGCATCCACGAGTTGTCTGTGGAGCA
 GCAGCTCTACTACAAGGAGATCACCGAGGCCTGCGTGGGCTCTGCGAGGCCAAGAGGGCGGAAGCCCTG
 CAAAGCATTGCCACGGACCCCTGGACTGTATCAGATGCTGCCACGGTTCAGTACCTTTATCTCGGAGGGGG
 TCCGTGTGAACGTGGTTCAGAACAACTGGCCCTACTCATCTACCTGATGCGTATGGTGAAGCGCTGAT
 GGACAACCCACGCTCTATCTAGAAAAATACGTCCATGAGCTGATTCCAGCTGTGATGACCTGCATCGTG
 AGCAGACAGTTGTGCCTGCGACCAGATGTGGACAATCACTGGGCACTCCGAGACTTTGCTGCCCGCTCG
 TGGCCAGATCTGCAAGCATTTTAGCACAACTAACAACATCCAGTCCCGGATCACCAGACCTTAC
 CAAGAGCTGGGTGGACGAGAAGACGCCCTGGACGACTCGTTATGGCTCCATCGCAGGCTTGCTGAGCTG
 GGACACGATGTTATCAAGACTCTGATTCTGCCCGGCTGCAGCAGGAAGGGGAGCGGATCCGCAGTGTGC
 TGGACGGCCCTGTGCTGAGCAACATTGACCGGATTGGAGCAGACCATGTGCAGAGCCTCCTGCTGAAACA
 CTGTGCTCCTGTTCTGGCAAAGCTGCGCCACCGCCTGACAATCAGGACGCCTATCGGGCAGAATTGCGG
 TCCCTTGGGCCCTCCTCTGCTCCAGGTGGTCAAGGCTCGGGCCAGGCTGCTCTGCAGGCTCAGCAGG
 TCAACAGGACCACTCTGACCATCACGACGCCCGGCCACGCTGACCCTCTCGCAGGCCCCACAGCCTGG
 CCCTCGCACCCCTGGCTTGCTGAAGTTCTGGTCCATCGCACTTCTGTCCAGACACTGGTGTCTGCA
 CGAGCGGCTGCCACCACAGCCTTCCCTCCTCAACCAAGTTATTGTAAATGTCATCGTCTCCAGCG
 CCCCATCCACCCAGCAGTCTGTCCCTCAGCACCTCGGCCCGGCTCAGGTTCCACCACCACTTCGCC
 CGTCAACCACCCGTCGCCAGCGTGCAGCCATCGTCAAGTTGGTCTCCACCGCCACCACCGCACCCCC
 AGCACTGCTCCCTCTGGTCTGGGAGTGTCCAGAAGTACATCGTGGTCTCACTTCCCCAACAGGGGAGG
 GCAAAGGAGGCCCCACCTCCCATCCTTCTCCAGTTCTCCCCGGCATCGTCCCCGTCCCCACTCAGCGG
 CAGTGCCCTTTGTGGGGGAAGCAGGAGGCTGGGACAGTCCCCCTCCAGCTCCAGGACTCCAAAAGCC
 AATGGCTCCAGCCCACTCCGGCTCCCTCAGCTGTCCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231390 representing NM_001190415
Red=Cloning site Green=Tags(s)


MTKRRLQTITACLQLREGSPSLHRGLHPSREEKRDSRMAEEKLKL SNTVLPSESMKVVAESMGIAQIQE
 ETCQLLTDEVSYRIKEIAQDALKFMMHGKRQKLTTSIDYALKLNVEPLYGFHAQEFIPFRFASGGGRE
 LYFYEEKEVDLSDIINTPLPRVPLDVCLKAHWLSIEGCQPAIPENPPPAPKEQQKAEATEPLKSAKPGQE
 EDGPLKGGKQGATTADGKGKEKKAPPLLEGAPLRKPRSIHEL SVEQQLYKYEITEACVGSCEAKRAEAL
 QSIATDPGLYQMLPRFSTFI SEGVRVNVVQNNLALLIYLMRMVKALMDNPTLYLEKYVHELIPAVMTCIV
 SRQLCLRPVDNHWALRDF AARLVAQICKHFSTTTNNIQSRITKFTKSWVDEKTPWTTRYGSIAGLAE
 GHDVIKTLILPRLQQEGERIRSVLDGPVLSNIDRIGADHVQSLLLKHCAPVLAKLRPPPDNDAYRAEFG
 SLGPLLCSQVVKARAQAALQAQVNRRTTLITQPRPTLSQAPQGPRTPLGLLKVPGSIALPVQTLVSA
 RAAAPPQSPPTKFI VMSSSSSAPSTQQVLSLSTSAPGSGSTTSPVTTTVPSVQPIVKLVSTATTAPP
 STAPSGPGSVQKIYVSLPPTGEGKGGPTSHSPVPPPASSPSPLSGSALCGGKQEAGDSPPPAPGTPKA
 NGSQPNSGSPQAP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001190415

ORF Size: 2142 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001190415.2
RefSeq ORF:	2145 bp
Locus ID:	6878
UniProt ID:	P49848
Cytogenetics:	7q22.1
Protein Families:	Transcription Factors
Protein Pathways:	Basal transcription factors
MW:	77.4 kDa
Gene Summary:	<p>Initiation of transcription by RNA polymerase II requires the activities of more than 70 polypeptides. The protein that coordinates these activities is transcription factor IID (TFIID), which binds to the core promoter to position the polymerase properly, serves as the scaffold for assembly of the remainder of the transcription complex, and acts as a channel for regulatory signals. TFIID is composed of the TATA-binding protein (TBP) and a group of evolutionarily conserved proteins known as TBP-associated factors or TAFs. TAFs may participate in basal transcription, serve as coactivators, function in promoter recognition or modify general transcription factors (GTFs) to facilitate complex assembly and transcription initiation. This gene encodes one of the smaller subunits of TFIID that binds weakly to TBP but strongly to TAF1, the largest subunit of TFIID. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2010]</p>

Product images:



Circular map for RC231390