

## Product datasheet for **RC231404**

### CD22 (NM\_001185099) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	CD22 (NM_001185099) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CD22
Synonyms:	SIGLEC-2; SIGLEC2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide  
Sequence:

>RC231404 representing NM\_001185099  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
GCCGCGATCGCC

ATGCATCTCCTCGGCCCTGGCTCCTGCTCCTGGTTCTAGAATACTTGGCTTCTCTGACTCAAGTAAAT  
GGGTTTTTGAGCACCTGAAACCCTCTACGCCTGGGAGGGGCCTGCGTCTGGATCCCCTGCACCTACAG  
AGCCCTAGATGGTGACCTGGAAGCTTTCATCCTGTTCCACAATCCTGAGTATAACAAGAACACCTCGAAG  
TTTGATGGGACAAGACTCTATGAAAGCACAAGGATGGGAAGGTTCTTCTGAGCAGAAAAGGTGCAAT  
TCCTGGGAGACAAGAATAAGAACTGCACACTGAGTATCCACCCGGTGCACCTCAATGACAGTGGTCACT  
GGGGCTGAGGATGGAGTCCAAGACTGAGAAATGGATGGAACGAATACACCTCAATGTCTCTGAAAGGCC  
TTTCCACCTCATATCCAGCTCCCTCCAGAAATCAAGAGTCCCAGGAAGTCACTCTGACCTGCTTGTCTGA  
ATTTCTCCTGCTATGGGTATCCGATCCAATTGCAGTGGCTCCTAGAGGGGTTCCAATGAGGAGGCTGC  
TGTCACCTCGACCTCCTTGACCATCAAGTCTGTCTTACCCTGGAGCGAGCTCAAGTTCTCCACAGTGG  
AGTCACCATGGGAAGATTGTGACCTGCCAGCTTCAAGGATGCAGATGGGAAGTTCCCTCCAATGACACGG  
TGCAGCTGAACGTGAAGCACACCCCGAAGTTGGAGATCAAGGTCACCTCCAGTGTGCCATAGTGAGGGA  
GGGGGACTCTGTGACCATGACCTGCGAGGTGAGCAGCAGCAACCCGGAGTACACGACGGTATCCTGGCTC  
AAGGATGGGACCTCGCTGAAGAAGCAGAATACATTCACGCTAAACCTGCGCGAAGTGACCAAGGACCAGA  
GTGGGAAGTACTGCTGTGAGGCTCCAATGACGTGGGCCCGGGAAGTTCGGAAGAAGTGTTCCTGCAAGT  
GCAGTATCCTCCCAAGAAGGTGACCACAGTATTCAAACCCATGCCGATTCGAGAAGGAGACACAGTGC  
ACCTTTTCTGTAACATAATCCAGTAAACCCAGTGTACCCGGTATGAATGAAACCCCATGGCCCTGCGC  
GGGAGGACCATCGCTTGGGGTCTGAAGATCCAAAACGTTGGCTGGGACAACACAACCATCGCCTGCGC  
AGCTTGTAAATAGTTGGTGTCTGTTGGGCTCCCTGTCGCTGAATGTCCAGTATGCCCCCGAGACGCTG  
AGGGTCCGGAATAAAGCCCTTCCGAGATTACTCTGGAACCTCGGTGAGCTCCAATGTGACTTCT  
CAAGCAGCACCCCAAGAAGTCCAGTCTTCTGGGAGAAAAATGGCAGGCTTCTGGGAAAGAAAGCCA  
GCTGAATTTTACTCCATCTCCCAGAAGATGCTGGGAGTTACAGCTGCTGGGTGAACAACCTCATAGGA  
CAGACAGCGTCCAAGGCTGGACACTTGAAGTGTGTATGCACCCAGGAGGCTGCGTGTGCTCATGAGCC  
CGGGGACCAAGTGTGAGGGGAAGAGTGCAACCCTGACCTGTGAGAGCGACGCCAACCTCCCGTCTC  
CCTACACCTGGTTTACTGGAATAACCAAAGCCTCCCTACCACAGCCAGAAGCTGAGATTGGAGCCG  
GTGAAGTCCAGCACTCGGTGCCTACTGGTGCAGGGGACCAACAGTGTGGGCAAGGGCGTTTCGCTC  
TCAGCACCTCACCGTCTACTATAGCCCGGAGACCATCGGCAGGCGAGTGGCTGTGGGACTCGGGTCTG  
CCTCGCCATCCTCATCTGGCAATCTGTGGGCTCAAGCTCCAGCGACGTTGGAAGAGGACACAGAGCCAG  
CAGGGGCTTACAGGAGAATTCCAGCGGCCAGAGCTTCTTTGTGAGGAATAAAAAGGTTAGAAGGGCCCCC  
TCTCTGAAGGCCCCACTCCCTGGGATGCTACAATCCAATGATGGAAGATGGCATTAGCTACACCACCT  
GCGCTTCCCGAGATGAACATACCACGAACTGGAGATGCAGAGTCTCAGAGATGCAGAGACCTCCCCG  
GACTGCGATGACACGGTCACTTATTCAGCATTGCACAAGCGCAAGTGGGCGACTATGAGAACGTCATTC  
CAGATTTCCAGAAGATGAGGGGATTACTACTCAGAGCTGATCCAGTTTGGGGTGGGGAGCGGCCCTCA  
GGCACAAGAAAAATGTGGACTATGTGATCCTCAAACAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTAA

**Protein Sequence:** >RC231404 representing NM\_001185099  
 Red=Cloning site Green=Tags(s)

MHLLGPWLLLLVLEYLAFSDSSKWVFEHPETLYAWEGACVWIPCTYRALDGDLESFILFHNPEYKNKNTSK  
 FDGTRL YESTKDGKVPSEQKRVQFLGDKNKNTLSIHPVHLNDSGQLGLRMEKTEKWMERIHLNVSERP  
 FPPHIQLPPEIQESQEVTLTCLLNFSYGYPIQLQWLLLEGVPMRQAAVTSTSLTIKSVFTRSELKFSQW  
 SHHGKIIVTCQLQDADGKFLSNDTVQLNVKHTPKLEIKVTPSDAIVREGDSVTMTCEVSSSNPEYTTVSWL  
 KDGTSLKKQNTFFLNLREVTKDQSGKYCCQVSNVDVGPGRSEEVFLQVQYPPKKVTTVIQNPMPPIREGDTV  
 TLSCNYSNSNPSVTRYEWKPHGAWEEP SLGVLKIQNVGWDNTTIIACAACNSWCSWASPVALNVQYAPRDV  
 RVRKIKPLSEIHSGNSVSLQCDFSSSHPKVEVQFFWEKNGRLLGKESQLNFDSISPEDAGSYSCWVNNSIG  
 QTASKAWTLEVL YAPRRLRVSMSPGDQVMEGKSATLTCESDANPPVSHYTFDWNQSLPYHSQKLRLEP  
 VKVQHSGAYWCQGTNSVGKGRSPLSTLVYYPETIGRRVAVGLGSCLA ILILAICGLKLRWRKRTQSQ  
 QGLQENSSGQSFVVRNKKVRRAPLSEGPLSLGCYNPMMEDGISYTTLRFPEMNIPRTGDAESSEMQRPP  
 DCDDTVTYSALHKRQVGDYENVIPDFPEDEGIHYSELIQFVGGERPQAQENVYVILKH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Chromatograms:** [https://cdn.origene.com/chromatograms/ja3495\\_c04.zip](https://cdn.origene.com/chromatograms/ja3495_c04.zip)

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF


**ACCN:** NM\_001185099

**ORF Size:** 2277 bp

<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_001185099.2</a>
<b>RefSeq Size:</b>	3036 bp
<b>RefSeq ORF:</b>	2280 bp
<b>Locus ID:</b>	933
<b>UniProt ID:</b>	<a href="#">P20273</a>
<b>Cytogenetics:</b>	19q13.12
<b>Protein Families:</b>	Druggable Genome, Transmembrane
<b>Protein Pathways:</b>	B cell receptor signaling pathway, Cell adhesion molecules (CAMs), Hematopoietic cell lineage
<b>MW:</b>	85.7 kDa

**Gene Summary:**

Mediates B-cell B-cell interactions. May be involved in the localization of B-cells in lymphoid tissues. Binds sialylated glycoproteins; one of which is CD45. Preferentially binds to alpha-2,6-linked sialic acid. The sialic acid recognition site can be masked by cis interactions with sialic acids on the same cell surface. Upon ligand induced tyrosine phosphorylation in the immune response seems to be involved in regulation of B-cell antigen receptor signaling. Plays a role in positive regulation through interaction with Src family tyrosine kinases and may also act as an inhibitory receptor by recruiting cytoplasmic phosphatases via their SH2 domains that block signal transduction through dephosphorylation of signaling molecules. [UniProtKB/Swiss-Prot Function]

**Product images:**


Circular map for RC231404