

Product datasheet for **RC231415LIV**

SATBI (NM_001195470) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	SATBI
Synonyms:	DEFDA; KTZSL
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001195470
ORF Size:	2385 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC231415).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_001195470.1](#)

RefSeq ORF: 2388 bp

Locus ID: 6304

UniProt ID: [Q01826](#)

Cytogenetics: 3p24.3

Protein Families: Druggable Genome, Transcription Factors



MW: 89.6 kDa

Gene Summary: This gene encodes a matrix protein which binds nuclear matrix and scaffold-associating DNAs through a unique nuclear architecture. The protein recruits chromatin-remodeling factors in order to regulate chromatin structure and gene expression. [provided by RefSeq, Apr 2016]