

## Product datasheet for **RC231429L4V**

### YY1AP1 (NM\_001198904) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	YY1AP1 (NM_001198904) Human Tagged ORF Clone Lentiviral Particle
Symbol:	YY1AP1
Synonyms:	GRNG; HCCA1; HCCA2; YY1AP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001198904
ORF Size:	2604 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231429).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001198904.1</a> , <a href="#">NP_001185833.1</a>
RefSeq ORF:	2607 bp
Locus ID:	55249
UniProt ID:	<a href="#">Q9H869</a>
Cytogenetics:	1q22
MW:	96 kDa
Gene Summary:	The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]



[View online »](#)