

Product datasheet for **RC231430L1V**

OXRI (NM_001198533) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Symbol:	OXRI
Synonyms:	CHEGDD; Nbla00307; TLDC3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001198533
ORF Size:	2619 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(RC231430).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq:	NM_001198533.1 , NP_001185462.1
RefSeq Size:	4716 bp
RefSeq ORF:	2622 bp
Locus ID:	55074
UniProt ID:	Q8N573
Cytogenetics:	8q23.1



MW: 97.9 kDa

Gene Summary: May be involved in protection from oxidative damage.[UniProtKB/Swiss-Prot Function]