

Product datasheet for **RC231434L4V**

MAP3K12 (NM_001193511) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | MAP3K12 (NM_001193511) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | MAP3K12 |
| Synonyms: | DLK; MEKK12; MUK; ZPK; ZPKP1 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_001193511 |
| ORF Size: | 2676 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC231434). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001193511.1 |
| RefSeq Size: | 3640 bp |
| RefSeq ORF: | 2679 bp |
| Locus ID: | 7786 |
| UniProt ID: | Q12852 |
| Cytogenetics: | 12q13.13 |
| Protein Families: | Druggable Genome, Protein Kinase |
| Protein Pathways: | MAPK signaling pathway |



[View online »](#)

MW: 96.8 kDa

Gene Summary: This gene encodes a member of the serine/threonine protein kinase family. This kinase contains a leucine-zipper domain and is predominately expressed in neuronal cells. The phosphorylation state of this kinase in synaptic terminals was shown to be regulated by membrane depolarization via calcineurin. This kinase forms heterodimers with leucine zipper containing transcription factors, such as cAMP responsive element binding protein (CREB) and MYC, and thus may play a regulatory role in PKA or retinoic acid induced neuronal differentiation. Alternatively spliced transcript variants encoding different proteins have been described.[provided by RefSeq, Jul 2010]