

Product datasheet for **RC231447L1V**

NMDAR1 (GRIN1) (NM_001185090) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | Lentiviral Particles |
| Product Name: | NMDAR1 (GRIN1) (NM_001185090) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | NMDAR1 |
| Synonyms: | GluN1; MRD8; NDHMSD; NDHMSR; NMD-R1; NMDA1; NMDAR1; NR1 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_001185090 |
| ORF Size: | 2829 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC231447). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_001185090.1 |
| RefSeq ORF: | 2832 bp |
| Locus ID: | 2902 |
| UniProt ID: | Q05586 |
| Cytogenetics: | 9q34.3 |
| Protein Families: | Druggable Genome, Ion Channels: Glutamate Receptors, Transmembrane |
| Protein Pathways: | Alzheimer's disease, Amyotrophic lateral sclerosis (ALS), Calcium signaling pathway, Huntington's disease, Long-term potentiation, Neuroactive ligand-receptor interaction |



[View online »](#)

MW: 106.5 kDa

Gene Summary: The protein encoded by this gene is a critical subunit of N-methyl-D-aspartate receptors, members of the glutamate receptor channel superfamily which are heteromeric protein complexes with multiple subunits arranged to form a ligand-gated ion channel. These subunits play a key role in the plasticity of synapses, which is believed to underlie memory and learning. Cell-specific factors are thought to control expression of different isoforms, possibly contributing to the functional diversity of the subunits. Alternatively spliced transcript variants have been described. [provided by RefSeq, Jul 2008]