

Product datasheet for **RC231478L3V**

WDR65 (CFAP57) (NM_001195831) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	WDR65 (CFAP57) (NM_001195831) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WDR65
Synonyms:	VWS2; WDR65
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001195831
ORF Size:	3849 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC231478).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001195831.2
RefSeq Size:	4180 bp
RefSeq ORF:	3852 bp
Locus ID:	149465
UniProt ID:	Q96MR6
Cytogenetics:	1p34.2
MW:	149 kDa


[View online »](#)

Gene Summary:

This protein encoded by this gene belongs to the WD repeat-containing family of proteins, which function in the formation of protein-protein complexes in a variety of biological pathways. This family member is thought to function in craniofacial development, possibly in the fusion of lip and palate. A missense mutation in this gene is associated with Van der Woude syndrome 2. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Aug 2011]