

Product datasheet for RC231506

TMEM127 (NM_001193304) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM127 (NM_001193304) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TMEM127
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC231506 representing NM_001193304 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTACGCCCCGGAGGCGCAGGGCTGCCGGCGGGCGCCGGCGGAGGAGCCGGGAGGCAGCGCTCTGC
CCAAGCAGCCGGAGCGTAGCCTGGCCTCGGCCCTGCCTGGCGCCCTGTCTATCACGGCGCTGTGCACTGC
CCTCGCCGAGCCCGCCTGGTTGCACATCCACGGAGGCACCTGTTTCGCGCCAGGAGCTGGGGTCTCCGAC
GTGTTGGGCTATGTGACCCGGACCTGCTGAAAGATTTCTGCATGAATCCCAGACAGTGTCTCTCTGC
GGGTATCGCCGCTTCTGTTTCTGGGCATCCTGTGTAGTCTCTCCGCTTTCCTTCTGGATGCTTTGG
GCCGAAGCATCCTGCTCTGAAGTCACTCGTCGCTATGCCTTCGCCCATATCCTAACGGTTCGAGTGT
GCCACCGTCATTGGCTTTTCTTATTGGGCTTCTGAACTCATCTTGGCCCAGCAGCAGCAGCATAAGAAGT
ACCATGGATCCCAGGTCTATGTCACCTTCGCCGTTAGCTTCTACCTGGTGGCAGGAGCTGGTGGAGCCTC
AATCCTGGCCACGGCAGCCAACCTCCTGCGCCACTACCCACAGAGGAAGAGGAGCAGGCCTGGAGCTG
CTCTCAGAGATGGAAGAGAACGAGCCCTACCCGGCGGAATATGAGGTCATCAACCAGTTCAGCCACCCC
CTGCTTACACCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC231506 representing NM_001193304
Red=Cloning site Green=Tags(s)

MYAPGGAGLPGGRRRRSPGGSALPKQPERSLASALPGALSITALCTALAEPAWLHIHGTCRSRQELGVSD
 VLGYYHPDLLKDFCMNPQTVLLLRVIAAFCLGILCSLSAFLLDVFQPKHPALKITRRYAFHILTVLQC
 ATVIGFSYWASELILAQQQHKYHGSQVYVTFVAVSFYLVAGAGGASILATAANLLRHYPTEEEEQALEL
 LSEMEENEPYPAEYEVINQFQPPPAYTP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/ja1479_c09.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001193304

ORF Size: 714 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001193304.3](#)

RefSeq Size: 4570 bp

RefSeq ORF: 717 bp

Locus ID: 55654

UniProt ID: [O75204](#)


Cytogenetics: 2q11.2

Protein Families: Transmembrane

MW: 25.8 kDa

Gene Summary: This gene encodes a transmembrane protein with 3 predicted transmembrane domains. The protein is associated with a subpopulation of vesicular organelles corresponding to early endosomal structures, with the Golgi, and with lysosomes, and may participate in protein trafficking between these structures. Mutations in this gene and several other genes cause pheochromocytomas. Alternatively spliced transcript variants encoding the same protein have been identified. [provided by RefSeq, Aug 2010]

Product images:



Circular map for RC231506