

Product datasheet for **RC231513**

C12orf62 (COX14) (NM_001257133) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	C12orf62
Synonyms:	C12orf62; MC4DN10; PCAG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC231513 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCCAACTGGCAAGCAGCTAGCTGACATTGGCTATAAGACCTTCTACCTCCATGATGCTTCTCACTG
TGTATGGGGGTACCTCTGCAGTGTCCGAGTCTACCACTATTTCCAGTGGCGCAGGGCCAGCGCCAGGC
CGCAGAAGAACAGAAGACCTCAGGAATCATG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231513 protein sequence
Red=Cloning site Green=Tags(s)

MPTGKQLADIGYKTFSTSMMLLTVYGGYLCVSRVYHYFQWRRARQAEEQKTSKIM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6412_d05.zip

Restriction Sites: SgfI-MluI



Cloning Scheme:



ACCN: NM_001257133

ORF Size: 171 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001257133.2
RefSeq Size:	814 bp
RefSeq ORF:	174 bp
Locus ID:	84987
UniProt ID:	Q96I36
Cytogenetics:	12q13.12
Protein Families:	Transmembrane
MW:	6.6 kDa
Gene Summary:	This gene encodes a small single-pass transmembrane protein that localizes to mitochondria. This protein may play a role in coordinating the early steps of cytochrome c oxidase (COX; also known as complex IV) subunit assembly and, in particular, the synthesis and assembly of the COX I subunit of the holoenzyme. Mutations in this gene have been associated with mitochondrial complex IV deficiency. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012]

Product images:



Circular map for RC231513