

Product datasheet for **RC231552A1V**

Human ZNF706 (NM_001267708) AAV Particle

Product data:

Product Type: AAV Particles

Product Name: Human ZNF706 (NM_001267708) AAV Particle

Tag: Myc-DDK

Symbol: ZNF706

Synonyms: HSPC038; PNAS-106; PNAS-113

Mammalian Cell Selection: None

Vector: pAAV-AC-Myc-DDK (PS100089)

ORF Nucleotide Sequence: >RC231552 representing NM_001267708
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTCGTGGACAGCAGAAAATTCAGTCTCAGCAGAAAAATGCCAAAAAGCAAGCTGGACAAAAGAAGA
AACAAAGGACATGACCAAAGGCTGCTGCCAAAGCTGCCTTAATATACCTGCACTGTCTGTAGGACACA
AATGCCAGACCCTAAGACCTTCAAGCAGCACTTTGAGAGCAAGCATCCTAAGACTCCACTTCTCCAGAA
TTAGCTGATGTTCAAGCA

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231552 representing NM_001267708
Red=Cloning site Green=Tags(s)

MARGQQKIQSQKNAKKQAGQKKQGHDKAAAKAALIYTCTVCRTQMPDPKTFKQHFESKHPKTPLPPE
LADVQA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human

Serotype: AAV-2

ACCN: NM_001267708

ORF Size: 228 bp

Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001267708.1</u> , <u>NP_001254637.1</u>
RefSeq Size:	2914 bp
RefSeq ORF:	231 bp
Locus ID:	51123
UniProt ID:	<u>Q9Y5V0</u>
Cytogenetics:	8q22.3
MW:	8.5 kDa