

Product datasheet for RC231566

Bcl rambo (BCL2L13) (NM_001270735) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bcl rambo (BCL2L13) (NM_001270735) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: BCL2L13
Synonyms: BCL-RAMBO; Bcl2-L-13; MIL1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC231566 representing NM_001270735
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGTCCTCTTCTACTGTGCCTCTGGGATTTCACTATGAAACAAAGTATGTTGTTCTCAGCTACTTGG
GACTCCTCTCTCAAGAGAAGCTGCAAGAGCAACATCTTTCCTCACCCCAAGGGTTCAACTAGATATAGC
TTCACAATCTCTGGATCAAGAAATTTATTTAAAAGTTAAAAGTAAATGAAGAAGAGCTAAAATCTCTG
GACAAAGAAATTTCTGAAGGGCACTGTGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231566 representing NM_001270735
Red=Cloning site Green=Tags(s)

MASSTVPLGFHYETKYVLSYLGLLSQEKLQEQLSSPQGVQLDIASQSLDQEILLKVKTEIEEELKSL
DKEISEGHCV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001270735

ORF Size: 240 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270735.1, NP_001257664.1</u>
RefSeq Size:	4761 bp
RefSeq ORF:	243 bp
Locus ID:	23786
Cytogenetics:	22q11.21
Protein Families:	Druggable Genome, Transmembrane
MW:	9.4 kDa
Gene Summary:	This gene encodes a mitochondrially-localized protein with conserved B-cell lymphoma 2 homology motifs. Overexpression of the encoded protein results in apoptosis. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012]