

Product datasheet for RC231637

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

UBE2V1 (NM 001257394) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: UBE2V1 (NM_001257394) Human Tagged ORF Clone

Tag: Myc-DDK
Symbol: UBE2V1

Synonyms: CIR1; CROC-1; CROC1; UBE2V; UEV-1; UEV1; UEV1A

Vector:pCMV6-Entry (PS100001)E. coli Selection:Kanamycin (25 ug/mL)

Cell Selection: Neomycin

ORF Nucleotide >RC231637 representing NM_001257394
Sequence: Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC

GCCGCGATCGCC

ATGAAGGAAGACTTGAACCTAGAAAATTTCACTGCAAAGACAATTTATGAAAAACCGAATATACAGCCTTA AAATAGAATGTGGACCTAAATACCCCAGAAGCACCCCCCTTTGTAAGATTTGTAACAAAAAATTAATATGAA TGGAGTAAATAGTTCTAATGGAGTGGTGGACCCAAGAGCCATATCAGTGCTAGCAAAATGGCAGAATTCA TATAGCATCAAAGTTGTCCTGCAAGAGCTTCGGCGCCCTAATGATGTCTAAAGAAAATATGAAACTCCCTC

AGCCGCCGAAGGACAGTGTTACAGCAAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT

ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231637 representing NM_001257394

Red=Cloning site Green=Tags(s)

MKEDLNLENFTAKTIYENRIYSLKIECGPKYPEAPPFVRFVTKINMNGVNSSNGVVDPRAISVLAKWONS

YSIKVVLQELRRLMMSKENMKLPQPPEGQCYSN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: Sgfl-Mlul






Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_001257394

ORF Size: 309 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info



UBE2V1 (NM_001257394) Human Tagged ORF Clone - RC231637

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001257394.2</u>

 RefSeq Size:
 2392 bp

 RefSeq ORF:
 312 bp

 Locus ID:
 7335

 UniProt ID:
 Q13404

 Cytogenetics:
 20q13.13

Protein Families: Druggable Genome, Transcription Factors

MW: 12.3 kDa

Gene Summary: Ubiquitin-conjugating E2 enzyme variant proteins constitute a distinct subfamily within the E2

protein family. They have sequence similarity to other ubiquitin-conjugating enzymes but lack the conserved cysteine residue that is critical for the catalytic activity of E2s. The protein encoded by this gene is located in the nucleus and can cause transcriptional activation of the human FOS proto-oncogene. It is thought to be involved in the control of differentiation by altering cell cycle behavior. Alternatively spliced transcript variants encoding multiple isoforms have been described for this gene, and multiple pseudogenes of this gene have been identified. Co-transcription of this gene and the neighboring upstream gene generates a

rare transcript (Kua-UEV), which encodes a fusion protein comprised of sequence sharing

identity with each individual gene product. [provided by RefSeq, Apr 2012]