

Product datasheet for RC231639

Epigen (EPGN) (NM_001270991) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Epigen (EPGN) (NM_001270991) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: Epigen
Synonyms: ALGV3072; EPG; PRO9904
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC231639 representing NM_001270991
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTTTGGGAGTTCCAATATCAGTCTATCTTTATTCAACGCAATGACAGCACTGACCGAAGAGGCAG
 CCGTGACTGTAACACCTCCAATCACAGCCAGCAAGCTGACAACATAGAAGGCCATAGCCTTGAAGTT
 CTCACACCTTTGCCTGGAAGATCATAACAGTTACTGCATCAACGGTGCTTGTGCATTCCACCATGAGCTA
 GAGAAAGCCATCTGCAGGTGTTTACTGGTTATACTGGAGAAAGGTGTCTAAAATTGAAATCGCCTTACA
 ATGTCTGTTCTGGAGAAAGACGACCACTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231639 representing NM_001270991
 Red=Cloning site Green=Tags(s)

MALGVPIISVYLLFNAMTALTEEAIVTVPPIITAQQADNIEGPIALKFSHLCLDHNSYINGACAFHHEL
 EKAICRCFTGYTGERCLKLSPYNNVCSGERRPL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001270991

ORF Size: 309 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001270991.2](#)

RefSeq Size: 2542 bp

RefSeq ORF: 312 bp

Locus ID: 255324

UniProt ID: [Q6UW88](#)


Cytogenetics: 4q13.3

Protein Families: Transmembrane

MW: 11.7 kDa

Gene Summary: The protein encoded by this gene is a member of the epidermal growth factor family. Members of this family are ligands for the epidermal growth factor receptor and play a role in cell survival, proliferation and migration. This protein has been reported to have high mitogenic activity but low affinity for its receptor. Expression of this transcript and protein have been reported in cancer specimens of the breast, bladder, and prostate. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2012]

Product images:



Circular map for RC231639