

Product datasheet for **RC231654**

p53 AIP1 (TP53AIP1) (NM_001251964) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: p53 AIP1 (TP53AIP1) (NM_001251964) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: TP53AIP1
Synonyms: P53AIP1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC231654 representing NM_001251964
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGATCTTCCTCTGAGGCGAGCTTCAGATCTGCTCAAGCTTCCTGCAGTGGGGCCAGGAGGCAGGGCC
TGGGCAGGGGAGACCAGAACCTCTCGGTGATGCCTCCGAATGGCAGGGCTCAGACACACACCTGGCTG
GGTAAGTCCCTGCAGTAAAACCGAGACGGTCTTTGCCTGCCACAGCCCGGGCAGACTCTGCTCTCAC
CGTGGTCCGACATCCCAAGTTTCAGACTCACCAGGACCCAGTGACAGCATCTGGGTCCTCAGAGCTGC
ATCGGGACTGTCCCAGTTCAGAGCATTGGACAGAGCTGGGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231654 representing NM_001251964
Red=Cloning site Green=Tags(s)

MGSSEASFRSAQASCSGARRQGLGRGDQNLVMPNGRAQHTHPGWVSPCSENDRGLLPATAPGRLLCSH
RGADIPSFQTHQDPVTASGSSELHADCPQFRALDRAGN

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001251964

ORF Size: 324 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001251964.1 , NP_001238893.1
RefSeq Size:	3158 bp
RefSeq ORF:	327 bp
Locus ID:	63970
UniProt ID:	Q9HCN2
Cytogenetics:	11q24.3
Protein Families:	Druggable Genome
Protein Pathways:	p53 signaling pathway
MW:	11.7 kDa
Gene Summary:	This gene is specifically expressed in the thymus, and encodes a protein that is localized to the mitochondrion. The expression of this gene is inducible by p53, and it is thought to play an important role in mediating p53-dependent apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011]