

Product datasheet for **RC231656**

C10orf57 (TMEM254) (NM_001270369) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C10orf57 (TMEM254) (NM_001270369) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: TMEM254
Synonyms: bA369J21.6; C10orf57
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC231656 representing NM_001270369
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGGGCGGGGACGGCTGGGGAGCTCCCAGCCAACTCTGTGTCTGGGTTTTTCAAGAGGATGGTGAAGT
CGGAGGCCAGCAATGCGGCTACCGCCTATTTCCGTCCAGCAAGACCCTTGCCCTCTCTGGTCACAGTCTT
GGGGCTGGGCTACTTCGCGTGGGTGTCTTCTGGCCTCAGAGTATCCCTTATCAGAACCTTGGGCCCTG
GGCCCTTCACTCAGTACTTGGTGGACCACCATCACCCCTCCTGTGCAATGGGTATTGGCTTGCCTGGC
TGATTCATGTGGGAGAGTCTTGTATGCCATAGTATTGTGCAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231656 representing NM_001270369
Red=Cloning site Green=Tags(s)

MGAGTAGELPANSVSGFFKRMVKSEASNAATAYFRPARPLPSLVTVLGLGYFAWVFWPQSIQNLGPL
GPFTQYLVDHHTLLCNGYWLAWLIHVGESLYAIVLCK

TRTRPLEQK**LISEEDLA**ANDILDYKDDDDK**V**

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001270369

ORF Size: 324 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001270369.1</u> , <u>NP_001257298.1</u>
RefSeq Size:	2152 bp
RefSeq ORF:	327 bp
Locus ID:	80195
UniProt ID:	<u>Q8TBM7</u>
Cytogenetics:	10q22.3
Protein Families:	Druggable Genome, Transmembrane
MW:	12.3 kDa