

Product datasheet for RC231659

Nurim (NRM) (NM_001270710) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nurim (NRM) (NM_001270710) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nurim
Synonyms:	NRM29
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC231659 representing NM_001270710 Red=Cloning site Blue=ORF Green=Tags(s) TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC ATGGCCCCGCACTGCTCCTGATCCCTGCTGCCCTCGCCTCTTTCATCCTGGCCTTTGGCACCGGAGTGGAGTTTCGTGCGCTTACCTCCCTTCGGCCACTTCTGGAGGGATCCCGGAGTCTGGTGGTCCGGGTATACTACCATGTGCTGGGGCTGGGCGAGCCTCTGGCCCTGAAGTCTCCCCGGGCTCTCAGACTCTTCTCCACCTGCGCCACCCAGTGTGTGTGGAGCTGCTGACAGTGTGTGGTGGTGCCTACCCTGGGCACGGACCGTCTCTCCTTGCTTTCCTCCTACCCTCTACCTGGGCCTGGCTCACGGCT ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC231659 representing NM_001270710 Red=Cloning site Green=Tags(s) MAPALLLIPAALASFILAFGTGVEFVRFTSLRPLLGGIPESGGPGILPCAGAGRASGPEVSPGSQTLPPAPPSVCGAADSAVGGAYPGHGSPPCFPYPLPGPSRA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Restriction Sites:	SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM_001270710

ORF Size: 327 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001270710.2](#)

RefSeq Size: 1196 bp

RefSeq ORF: 330 bp
 Locus ID: 11270
 UniProt ID: [Q8IXM6](#)
 Cytogenetics: 6p21.33
 Protein Families: Transmembrane
 MW: 11 kDa

Gene Summary: The protein encoded by this gene contains transmembrane domains and resides within the inner nuclear membrane, where it is tightly associated with the nucleus. This protein shares homology with isoprenylcysteine carboxymethyltransferase enzymes. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Jul 2012]

Product images:



Circular map for RC231659