

## Product datasheet for RC231665A1V

### Human TRPC5OS (NM\_001195577) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human TRPC5OS (NM\_001195577) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** TRPC5OS  
**Synonyms:** TRPC5-AS1  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC231665 representing NM\_001195577  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGATTCTGTGTTAATTCATGTA CTACTCATTGATGGACTTGTGCTTGTGTAGCCAGTTAATAAGAATAG  
 CTGATGAGCTTTTACAATTCATTCTACAAGTACAAGAAGTTCCTTATGTAGAAGAAAATGGTAGAGCAGA  
 AGAGACTGAAGCAGATGCACCTTTCCCGAGGAGCCTTCGCTACCTGATCTCCCTGATCTCTCAGACTTA  
 GACTCAATACTTACACCAAGAGAGGATGAAGACCTAATATTTGATATAGATCAGGCTATGTTAGACATGG  
 ATAACCTATATGAAGATACAGTCTCTGGTATAAATGATGACTTAACAGGTGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC231665 representing NM\_001195577  
 Red=Cloning site Green=Tags(s)

MDSVLIHVLIDGLVACVAQLIRIADELLQFILQVQVEVPYVEENGRAEETEADAPLPEEPSLPDL.PDL.SDL  
 DSILTPREDEDLIFDIDQAMLMDNLYEDTVSGINDDLTDG

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001195577  
**ORF Size:** 333 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u>NM_001195577.1, NP_001182506.1</u>
<b>RefSeq Size:</b>	2816 bp
<b>RefSeq ORF:</b>	336 bp
<b>Locus ID:</b>	100329135
<b>UniProt ID:</b>	<u>A6NMA1</u>
<b>Cytogenetics:</b>	Xq23
<b>MW:</b>	12.3 kDa