

Product datasheet for RC231693


Calcipressin 3 (RCAN3) (NM_001251984) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Calcipressin 3 (RCAN3) (NM_001251984) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Calcipressin 3
Synonyms:	DSCR1L2; hRCN3; MCIP3; RCN3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC231693 representing NM_001251984 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC ATGTCCGGCGAAGTGC GG GACAAGTCCTATCTCCTGCCGCCCCAGCCTGTCAAGCAGTTCCTCATCTCCCCTCCAGCCTCTCCCCGGTGGGGTGAAGCAGAGCGAAGATGCGATGCCTGTTATAAATTATGATTTACTCTGTGCTGTTTCAAATTGGGACCAGGAGAGAAATATGAACCTCACGCGGGAACAGAGTCGACACCCAGCGTGGTGGTTCATGTCTGTGAAAGTGAACTGAAGAGGAAGAAGAGACAAAAACCCAAACAGAAAATTGCCAGACGAGGCGCCCGACCCTCCGACCGCAGCGTTGAATGAGCCCCAGACCTTTGATTGCGCGCTG ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC231693 representing NM_001251984 Red=Cloning site Green=Tags(s) MSGEVRDKSYLLPPQPVKQFLISPPASPPVGWKQSEDAMPVINYDLLCAVSKLGPGEKYELHAGTESTPSVVVHVCESETEEEEEETKNPKQKIAQTRRPDPPTAALNEPQTFDCAL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Restriction Sites:	SgfI-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM_001251984

ORF Size: 348 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001251984.1](#), [NP_001238913.1](#)

RefSeq Size: 2392 bp

RefSeq ORF:	351 bp
Locus ID:	11123
UniProt ID:	Q9UKA8
Cytogenetics:	1p36.11
MW:	13.2 kDa
Gene Summary:	Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC231693