

Product datasheet for **RC231704**

ZNF248 (NM_001267605) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZNF248 (NM_001267605) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ZNF248
Synonyms: bA162G10.3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC231704 representing NM_001267605
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAACAATCCCAGGAACAAGTGTCAAGGATGTATGTGTGGACTTCACTCAGGAAGAGTGGTATC
TGCTGGACCCTGCTCAGAAGATTCTATACAGAGATGTGATCCTGGAAAATTATAGCAATCTTGTCTCAGT
AGGGTATTGCATTACTAAACCAGAAGTGATCTTTAAGATCGAGCAAGGAGAAGAGCCCTGGATATTAGAA
AAAGGATTCCTCAAGCCAGTGCCACCCAGAAAGGAAATGGAAAGTTGATGACGTGTAGAGAGCAGCCAGG
AAAATGAAGATGACCATTTTGGGAGCTTCTATCCACAACAACAAAACATGGAATCTGAAGCTTATAAA
AGTTGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231704 representing NM_001267605
Red=Cloning site Green=Tags(s)
MNKSQEQVSFKDVCVDFTEEWYLLDPAQKILYRDVILENYSNLVSVGVCITKPEVIFKIEQGEEPWILE
KGFPSQCHPERKWKVDDVLESSQENEDDHFWEELFHNNKTWNLKLIKVE


TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV**

Restriction Sites: Sgfl-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:

ACCN:

NM_001267605

ORF Size:

357 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001267605.1 , NP_001254534.1
RefSeq Size:	3479 bp
RefSeq ORF:	360 bp
Locus ID:	57209
UniProt ID:	Q8NDW4
Cytogenetics:	10p11.21
Protein Families:	Transcription Factors
MW:	14.7 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]