

Product datasheet for RC231736

TRAPPC6A (NM_001270892) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRAPPC6A (NM_001270892) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	TRAPPC6A
Synonyms:	TRS33
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC231736 representing NM_001270892 Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC ATGGCGGATACTGTGTTGTTGAGTTTCTTCACACGGAGATGGTGGCTGAGCTGTGGGCTCACGACCCCGACCCCGCGCCGGGGTGAGCGCCGGGCTCCGTGGGGAGGAAGCGGGGCCACCAAGGCTGCCCCGGGAGACGCTGGCCTTCAGGGAGGAGCTGGATGTCCTCAAGTTCTTGCAAGACCTGTGGGTGGCGGTGTTCCAGAAGCAGATGGACAGCCTGCGACCAATCACCAGGGGACCTACGTCCTGCAAGACAACAGCTTCCCCCTCTCCTCCCGATGGCCTCTGGCCTGCAGTATCTGGAGGAAGCACCAAGTTCCTGGCCTTCACCTGCGGCCCTCTGCGCGGCCCTCTATACCCTGGGCAT ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
Protein Sequence:	>RC231736 representing NM_001270892 Red=Cloning site Green=Tags(s) MADTVLFEFLHTEMVAELWAHDPDPGPVSAGLRGEEAGATKAAPGDAGLQGGAGCPQVLVQRPVGGGVPEADGQPAHQSPGDLRPARQQLPPPPDGLWPAVSGGSTQVPLHLRPPARRPLYPGH TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Restriction Sites:	Sgfl-MluI


[View online »](#)

Cloning Scheme:



ACCN: NM 001270892

ORF Size: 381 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_001270892.2](#)

RefSeq Size: 770 bp

RefSeq ORF: 384 bp
 Locus ID: 79090
 UniProt ID: [O75865](#)
 Cytogenetics: 19q13.32
 MW: 13.4 kDa

Gene Summary: This gene encodes a component of the trafficking protein particle complex, which tethers transport vesicles to the cis-Golgi membrane. Loss of expression of the related gene in mouse affects coat and eye pigmentation, suggesting that the encoded protein may be involved in melanosome biogenesis. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2012]

Product images:



Circular map for RC231736