

## Product datasheet for **RC231743A1V**

### Human PFD6 (PFDN6) (NM\_001265595) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human PFD6 (PFDN6) (NM\_001265595) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** PFD6  
**Synonyms:** H2-KE2; HKE2; KE-2; PFD6  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC231743 representing NM\_001265595  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGCGGAGCTGATCCAGAAGAAGCTACAGGGAGAAGTGGAGAAATATCAACAGCTACAGAAGGACTTAA  
GTAATCCATGTCGGGGAGGCAGAACTTGAAGCACAATAACAGAAAATAATATCGTGAAAGAGGAACT  
GGCCCTGCTGGATGGTCCAACGTGGTCTTTAACTTCTGGTCCGGTCTAGTCAAACAGGAGCTGGGG  
GAGGCTCGGGCCACAGTAGGGAAGAGGCTGGACTATATCACAGCTGAAATTAAGCGATACGAATCCCAGC  
TTCGGGATCTTGAGCGGCAGTCAGAGCAACAGAGGGAGACCCCTTGCTCAGCTGCAGCAGGAGTCCAGCG  
GCCCCAGGCAGCAAAGGCAGGGGCTCCTGGCAAGGCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC231743 representing NM\_001265595  
Red=Cloning site Green=Tags(s)

MAELIQKQLQGEVEKYQQLQKDLKSKSMSGRQKLEAQLTENNIVKEELALLDGSNVVFKLLGPVLVKQELG  
EARATVGKRLDYITAEIKRYESQLRDLERQSEQQRETLAQLQEFQRAQAAKAGAPGKA

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_001265595



[View online »](#)

<b>ORF Size:</b>	387 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_001265595.1</a></u>
<b>RefSeq Size:</b>	597 bp
<b>RefSeq ORF:</b>	390 bp
<b>Locus ID:</b>	10471
<b>UniProt ID:</b>	<u><a href="#">O15212</a></u>
<b>Cytogenetics:</b>	6p21.32
<b>MW:</b>	14.6 kDa