

Product datasheet for RC231756

C18orf21 (NM_001201474) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: C18orf21 (NM_001201474) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: C18orf21
Synonyms: HsT3108; PNAS-124; PNAS-131; XTP13
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC231756 representing NM_001201474
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGC

ATGGTGAAAAAGTTCAAAGACTCCAAAAGTGTATTGTTGATCACTTGTAACATGCAACAGAACAGTGA
 AACATCATGGTAAAAGTAGAAGCTTTGTGTCAACATTGAAGAGCAATCCTGCCACTCCTACAAGTAACT
 CAGCCTGAAGACACCAGAGAGAAGGACTGCAACCCAAATCATGACATGTCTGGCTCGAAAGGCAAGAGC
 CCAGCATCGGTTTTCAGAACACCTACATCTGGACAGTCAGTATCTACTTGCTCCTCAAAGAACACCAGCA
 AAACAAAGAAACACTTCTCTCACTAAAAATGTTACTTAGTCAGAATGAATCCCAAAGATTCAAAGGT
 GGACTTCAGAAATTTCTTATCTTCTCTGAAGGGTGGACTTTTAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231756 representing NM_001201474
Red=Cloning site Green=Tags(s)

MVKKFKDSKSVLLITCKTCNRTVKHHGKSRSFVSTLKSNPATPTSLSLKTERRTANPNHDMSGSKGKS
 PASVFRPTSGQSVSTCSSKNTSKTKKHFSQLKMLLSQNESQKIPKVDFRNLSSLKGGLLK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


[View online »](#)

Cloning Scheme:


ACCN: NM_001201474

ORF Size: 396 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001201474.1](#), [NP_001188403.1](#)

RefSeq Size: 1069 bp

RefSeq ORF: 399 bp
 Locus ID: 83608
 UniProt ID: Q32NC0
 Cytogenetics: 18q12.2
 MW: 15 kDa

Product images:



Circular map for RC231756