

## Product datasheet for RC231757

### C18orf21 (NM\_001201475) Human Tagged ORF Clone

#### Product data:

|                           |                                                                                |
|---------------------------|--------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                            |
| Product Name:             | C18orf21 (NM_001201475) Human Tagged ORF Clone                                 |
| Tag:                      | Myc-DDK                                                                        |
| Symbol:                   | C18orf21                                                                       |
| Synonyms:                 | HsT3108; PNAS-124; PNAS-131; XTP13                                             |
| Mammalian Cell Selection: | Neomycin                                                                       |
| Vector:                   | pCMV6-Entry (PS100001)                                                         |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                           |
| ORF Nucleotide Sequence:  | >RC231757 representing NM_001201475<br>Red=Cloning site Blue=ORF Green=Tags(s) |

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGGTGAAAAAGTTCAAAGACTCCAAAAGTGTATTGTTGATCACTTGTAACATGCAACAGAACAGTGA  
 AACATCATGGTAAAAGTAGAAGCTTTGTGTCAACATTGAAGAGCAATCCTGCCACTCCTACAAGTAACT  
 CAGCCTGAAGACACCAGAGAGAAGGACTGCAACCCAAATCATGACATGTCTGGCTCGAAAGGCAAGAGC  
 CCAGCATCGGTTTTCAGAACACCTACATCTGGACAGTCAGTATCTACTTGCTCCTCAAAGAACACCAGCA  
 AAACAAAGAAACACTTCTCTCACTAAAAATGTTACTTAGTCAGAATGAATCCCAAAGATTCAAAGGT  
 GGACTTCAGAAATTTCTTATCTTCTCTGAAGGGTGGACTTTTAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

|                   |                                                                       |
|-------------------|-----------------------------------------------------------------------|
| Protein Sequence: | >RC231757 representing NM_001201475<br>Red=Cloning site Green=Tags(s) |
|-------------------|-----------------------------------------------------------------------|

MVKKFKDSKSVLLITCKTCNRTVKHHGKSRSFVSTLKSNPATPTSKLSLKTPEPRTANPNHDMSGSKGKS  
 PASVFRPTSGQSVSTCSSKNTSKTKKHFSQLKMLLSQNESQKIPKVDFRNLSSLKGGLLK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

|                    |           |
|--------------------|-----------|
| Restriction Sites: | SgfI-MluI |
|--------------------|-----------|


[View online »](#)

**Cloning Scheme:**


**ACCN:** NM\_001201475

**ORF Size:** 396 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.


**Note:** Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:** [NM\\_001201475.1](#), [NP\\_001188404.1](#)

**RefSeq Size:** 1136 bp

RefSeq ORF: 399 bp  
 Locus ID: 83608  
 UniProt ID: [Q32NC0](#)  
 Cytogenetics: 18q12.2  
 MW: 15 kDa

## Product images:



Circular map for RC231757