

Product datasheet for **RC231824**

RPL17 (NM_001199345) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPL17 (NM_001199345) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: RPL17
Synonyms: L17; PD-1; RPL23
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC231824 representing NM_001199345
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCATATACGAAAAGCCACGAAGTATCTGAAAGATGTCACCTTACAGAAACAGTGTGTACCATTCCGAC
GTTACAATGGTGGAGTTGGCAGGTGTGCGCAGGCCAAGCAATGGGGCTGGACACAAGGTCGGTGGCCAA
AAAGAGTCTGAATTTTTGCTGCACATGCTTAAAAACGCAGAGAGTAATGCTGAACCTAAGGTTTAGAT
GTAGATTCTCTGGTCATTGAGCATATCCAAGTGAACAAAGCACCTAAGATGCGCCGCCGACCTACAGAG
CTCATGGTCGGATTAACCCATACATGAGCTCTCCCTGCCACATTGAGATGATCCTTACGGAAAAGGAACA
GATTGTTCTAAACCAGAAGAGGAGTTGCCAGAAGAAAAGATATCCCAGAAGAACTGAAGAAACAA
AACTTATGGCACGGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC231824 representing NM_001199345
Red=Cloning site Green=Tags(s)
MHIRKATKYLKDVTLQKQCVFRRYNGGVGRCAQAKQWGWQGRWPKSAEFLHMLKNAESNAELKGLD
VDSLVIIEHIQVNKAPKMRRTYRAHGRINPYMSSPCHIEMILTEKEQIVPKPEEEVAQKKKISQKLLKKQ
KLMARE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001199345

ORF Size: 438 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001199345.2
RefSeq Size:	695 bp
RefSeq ORF:	441 bp
Locus ID:	6139
UniProt ID:	P18621
Cytogenetics:	18q21.1
Protein Pathways:	Ribosome
MW:	17.5 kDa
Gene Summary:	Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 60S subunit. The protein belongs to the L22P family of ribosomal proteins. It is located in the cytoplasm. This gene has been referred to as rpl23 because the encoded protein shares amino acid identity with ribosomal protein L23 from <i>Halobacterium marismortui</i> ; however, its official symbol is RPL17. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring downstream C18orf32 (chromosome 18 open reading frame 32) gene. [provided by RefSeq, Dec 2010]