

Product datasheet for RC231943

PBR (TSPO) (NM_001256531) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PBR (TSPO) (NM_001256531) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PBR
Synonyms: BPBS; BZRP; DBI; IBP; MBR; mDRC; PBR; PBS; pk18; PKBS; PTBR
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC231943 representing NM_001256531.
Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCCCGCCCTGGGTGCCCGCATGGGCTTACGCTGGCGCCAGCCTGGGGTCTTCGTGGGCTCC
CGCTTTGTCCACGGCGAGGGTCTCCGCTGGTACGCCGGCCTGCAGAAGCCCTCGTGGCACCCGCCAC
TGGGTGCTGGGCCTGTCTGGGGCAGCTCTACTAGCCATGGGGTACGGCTCTACCTGGTCTGGAAA
GAGCTGGGAGGCTTACAGAGAAGGCTGTGGTCCCTGGGCCTCTACACTGGGCAGCTGGCCCTGAAC
TGGGATGGCCCCCATCTTCTTTGGTGCCGACAAATGGGCTGGGCCTGGTGGATCTCCTGCTGGTC
AGTGGGGCGGCGCAGCCACTACCGTGGCCTGGTACCAGGTGAGCCCGCTGGCCGCCCGCTGCTCTAC
CCCTACCTGGCCTGGCTGGCCTTCGCGACCACACTCAACTACTGCGTATGGCGGGACAACCATGGCTGG
CATGGGGGACGGCGGCTGCCAGAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Protein Sequence: >Peptide sequence encoded by RC231943
Blue=ORF Red=Cloning site Green=Tag(s)

```
MAPPWVPAMGFTLAPSLGCFVGSRFVHGEGLRWYAGLQKPSWHPHVLGPVWGTLYSAMGYGSYLWVK
ELGGFTEKAVVPLGLYTGQLALNWAAPPFFGARQMGWALVDLLLVSAAAATTVAWYQVSPLAARLLY
PYLAWLAFATTLNYCVWRDNHGWHGRRLPE
TRTRPLEQKLISEEDLAANDILDYKDDDDKV
```

Chromatograms: https://cdn.origene.com/chromatograms/mk6385_d06.zip

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_001256531

ORF Size: 507 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq Size: 1025 bp

RefSeq ORF: 510 bp

Locus ID: 706


UniProt ID: [P30536](#)

Cytogenetics: 22q13.2

Protein Families: Druggable Genome, Transmembrane
Protein Pathways: Neuroactive ligand-receptor interaction
MW: 18.8 kDa

Gene Summary: Present mainly in the mitochondrial compartment of peripheral tissues, the protein encoded by this gene interacts with some benzodiazepines and has different affinities than its endogenous counterpart. The protein is a key factor in the flow of cholesterol into mitochondria to permit the initiation of steroid hormone synthesis. Alternatively spliced transcript variants have been reported; one of the variants lacks an internal exon and is considered non-coding, and the other variants encode the same protein. [provided by RefSeq, Feb 2012]

Product images:



Circular map for RC231943