

Product datasheet for **RC232016**

COPZ1 (NM_001271736) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: COPZ1 (NM_001271736) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: COPZ1
Synonyms: CGI-120; COPZ; HSPC181; zeta-COP; zeta1-COP
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC232016 representing NM_001271736
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGAGGCTTCAGGACTGAAGGTATGTTCTGATCTCTTCAGGAACCTCCCTGTATACTGTCAAAGCCA
TCCTGATTCTGGACAATGATGGAGATCGACTTTTTGCCAAGTACTATGACGACACCTACCCAGTGTCAA
GGAGCAAAGGCCTTTGAGAAGAACATTTCAACAAGACCCATCGGACTGACAGTGAAATTGCCCTTTG
GAAGGCCTGACAGTGGTATACAAAAGCAGTATAGATCTCTATTTCTATGTGATTGGCAGCTCCTATGAAA
ATGAGCTGATGCTTATGGCTGTTCTGAACTGTCTCTTCGACTCATTGAGCCAGATGCTGAGGAAAAATGT
AGAAAAGCGAGCACTGCTGGAGAACATGGAGGGCTGTTCTTGGCTGTGGATGAAATTGTAGATGGAGGG
GTGATCCTAGAGAGTGATCCCCAGCAGGTGACACCGGGTGGCATTAAAGGGTGAAGATGTCCCCCTTA
CGGAGCAGACCGTGTCTCAGGTGCTGCAGTCAGCCAAAGAACAGATCAAGTGGTCACTCCTTCGG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC232016 representing NM_001271736
Red=Cloning site Green=Tags(s)

MGGFRTEGMFVSLQEPSLYTVKAILILDNDGDRLFACYDDTYPVKEQKAFKFNIFNKTHRTDSEIALL
EGLTVVYKSSIDLYFYVIGSSYENELMLMAVLNCLFDSLSQLRKNVEKRALLENMEGLFLAVDEIVDGG
VILESDPQQVVRVALRGEDVPLTEQTVSQVLQSAKEQIKWLLR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI



Cloning Scheme:


ACCN: NM_001271736

ORF Size: 555 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271736.2](#)

RefSeq Size: 2038 bp

RefSeq ORF: 558 bp

Locus ID: 22818


UniProt ID: [P61923](#)

Cytogenetics: 12q13.13

MW: 21.5 kDa

Gene Summary: This gene encodes a subunit of the cytoplasmic coatamer protein complex, which is involved in autophagy and intracellular protein trafficking. The coatamer protein complex is comprised of seven subunits and functions as the coat protein of coat protein complex (COP)I-vesicles. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Nov 2012]

Product images:



Circular map for RC232016