

Product datasheet for RC232040

SLC39A9 (NM_001252151) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: SLC39A9 (NM_001252151) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: SLC39A9
Synonyms: ZIP-9; ZIP9
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232040 representing NM_001252151
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTTGCTGGTGGACCAGATTGGTAACTCCCATGTGCATTCTACTGACGATCCAGAAGCAGCAAGGTCTA
GCAATTCAAAATCACCACCACGCTGGGTCTGGTTGCCATGCTGCAGCTGATGGTGTGCTTTGGGAGC
AGCAGCATCTACTTCACAGACCAGTGTCCAGTTAATTGTGTTTGTGGCAATCATGCTACATAAGGCACCA
GCTGCTTTTGGACTGGTTTCCTTCTTGATGCATGCTGGCTTAGAGCGGAATCGAATCAGAAAGCACTTGC
TGGTCTTTGCATTGGCAGCACCAGTTATGTCCATGGTGACATACTTAGGACTGAGTAAGAGCAGTAAAGA
AGCCCTTTCAGAGGTGAACGCCACGGGAGTGGCCATGCTTTTCTCTGCCGGGACATTTCTTTATGTTGCC
ACAGTACATGTCCCTCCCTGAGGTGGGCGGAATAGGGCACAGCCACAAGCCCGATGCCACGGGAGGGAGAG
GCCTCAGCCGCCTGGAAGTGGCAGCCCTGGTCTGGGTTGCCTCATCCCTCATCCTGTCTAGTAGGACA
CCAGCAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC232040 representing NM_001252151
Red=Cloning site Green=Tags(s)

MLLVLDQIGNSHVHSTDDPEAARSSNSKITTTTLGLVVHAAADGVALGAAASTSQTSVQLIVFVAIMLHKAP
AAFGLVSFLMHAGLERNRIRKHLVLFALAAPVMSMVTYLGLSKSKEALSEVNATGVAMLF SAGTFLYVA
TVHVLPEVGGIGHSHKPDATGGRGLSRLEVAALVLGCL IPLILSVGHQH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001252151

ORF Size: 567 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001252151.2](#)

RefSeq Size: 4865 bp

RefSeq ORF: 570 bp

Locus ID: 55334

UniProt ID: [Q9NUM3](#)

Cytogenetics: 14q24.1

Protein Families: Transmembrane

MW: 20 kDa

Gene Summary: May act as a zinc-influx transporter.[UniProtKB/Swiss-Prot Function]