

## Product datasheet for **RC232076**

### SCNM1 (NM\_001204856) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** SCNM1 (NM\_001204856) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** SCNM1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC232076 representing NM\_001204856  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCTTCGGGATGGACGCTTTGCTTGTCATCTGCCCATCGACCCGACTGGACACCCTGGCCATGC  
TGACTGCCACCCTGCAGGCAAGAAACATCTGTCCAGCTTGACAGCTTTTCTATGGCAAGAAGCAGCCGG  
AAAGGAAAGAAAGCAGAATCCAAAACATCAGAATGAATTGAGAAGGGAAGAAACAAAGCTGAGGCTCCT  
CTGCTAACTCAGACACGACTTATCACCCAGAGTGCTCTGCACAGAGCTCCCACTATAACAGTTGCTGCC  
GCCGGAAGTACAGACCAGAAGCCCTGGTCCCTCTGTCTCCCTTTCCCTATGCCACCCTCAGAGGTCAA  
ACTCAAAGTGGGAAGATCAGTAGGGAACCTGAACCTGCGGCTGGCCACAGGCCGAGGAGTCAGCAACT  
GTCTCAGCCCTGCACCCATGAGCCCAAGAAGACGAGCCCTGGACCATTATCTACCCTTCGAAGCT  
CTGGATGGATCCCAGATGGACGAGGTCGATGGTAAAAGATGAAAATGTTGAGTTTGACTCTGATGAGGA  
GGAACCACCTGATCTCCCTTGGAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC232076 representing NM\_001204856  
Red=Cloning site Green=Tags(s)

MLRDGRFACAIICPHRPVLDLAMLTAHRAGKKHLSSLQLFYGKKQPGKERKQNPKHQNELRREETKAEAP  
LLTQTRLITQSALHRAPHYNSCCRRKYRPEAPGPSVSLSPMPPSEVKLQSGKISREPEPAAGPQAEESAT  
VSAPAPMSPTRRRALDHYLTLRSSGWIPDGRGRWVKDENVEFDSDEEPPDLPLD


**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001204856

ORF Size: 585 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001204856.2</a></u>
<b>RefSeq Size:</b>	2085 bp
<b>RefSeq ORF:</b>	588 bp
<b>Locus ID:</b>	79005
<b>UniProt ID:</b>	<u><a href="#">Q9BWG6</a></u>
<b>Cytogenetics:</b>	1q21.3
<b>MW:</b>	22.3 kDa
<b>Gene Summary:</b>	SCNM1 is a zinc finger protein and putative splicing factor. In mice, Scnm1 modifies phenotypic expression of Scn8a (MIM 600702) mutations (Buchner et al., 2003 [PubMed 12920299]).[supplied by OMIM, Oct 2009]