

## Product datasheet for **RC232092**

### IL18BP (NM\_005699) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** IL18BP (NM\_005699) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** IL18BP  
**Synonyms:** FVH; IL18BPa  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >RC232092 representing NM\_005699  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGACCATGAGACACAAC**TGGACACCAGACCTCAGCC**TTTGTGGGTCTGCTCTGTGTGCCACGTCTG  
TCACTCTCTGGTCAGAGCCACACCTGTCTCGCAGACCACCACAGCTGCCACTGCCTCAGTTAGAAGCAC  
AAAGGACCCCTGCCCTCCAGCCCCAGTGTCCAGCAGCTAAGCAGTGTCCAGCATTGGAAGTGACC  
TGGCCAGAGGTGGAAGTGCCACTGAATGGAACGCTGAGCTTATCCTGTGTGGCCTGCAGCCGTTCCCA  
ACTTCAGCATCCTCTACTGGCTGGCAATGGTTCCTTCATTGAGCACCTCCAGGCCGACTGTGGGAGGG  
GAGCACCAGCCGGAACGTGGGAGCACAGGTACGACGCTGTGCAAGGCCTTGGTGTGGAGCAGCTGACC  
CCTGCCCTGCACAGCACTTCTCCTGTGTGCTCGTGGACCCTGAACAGGTTGTCCAGCGTCACGTCTG  
TCTGGCCAGCTCTGGGTGAGGAGCCCAAGGAGAGGCCTCCAGGAACAGGAGGACTCTGCTTCCATAT  
GTGGGGAGGAAAGGGTGGGCTCTGCCAGAGCAGCCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC232092 representing NM\_005699  
Red=Cloning site Green=Tags(s)

MTMRHNWTPDLSPLWVLLLCAHVVTLLVRATPVSQTTTAAATASVRSTKDPSPQPPVFAAKQCPALEVT  
WPEVEVPLNGTSLSCVACSRFPNFSILYWLNGSFIHLPGRLEWEGSTSRERSTGTQLCKALVLEQLT  
PALHSTNFSCVLVDPEQVVQRHVLAQLWVRSRRGLQEQEELCFHMWGGKGLCQSSL

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM\_005699

ORF Size: 597 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
  2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
  3. Close the tube and incubate for 10 minutes at room temperature.
  4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
  5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_005699.3](#), [NP\\_005690.2](#)

RefSeq Size: 2796 bp

RefSeq ORF: 600 bp


Locus ID: 10068

UniProt ID: [O95998](#)

**Cytogenetics:** 11q13.4  
**Protein Families:** Druggable Genome, Secreted Protein  
**MW:** 22.4 kDa

**Gene Summary:** The protein encoded by this gene functions as an inhibitor of the proinflammatory cytokine, IL18. It binds IL18, prevents the binding of IL18 to its receptor, and thus inhibits IL18-induced IFN-gamma production, resulting in reduced T-helper type 1 immune responses. This protein is constitutively expressed and secreted in mononuclear cells. Elevated level of this protein is detected in the intestinal tissues of patients with Crohn's disease. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Feb 2011]

**Product images:**



Circular map for RC232092