

## Product datasheet for **RC232110**

### YIPF5 (NM\_001271732) Human Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** YIPF5 (NM\_001271732) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** YIPF5  
**Synonyms:** FinGER5; SB140; SMAP-5; SMAP5; YIP1A  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC232110 representing NM\_001271732  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGATGCAGCCACAACAGCCATACACCGGGCAGATTTACCAGCCAACCTCAGGCATATACTCCAGCTTCAC  
CTCAGCCTTTCTATGAAACAACCTTTGAGGATGAGCCACCTTTATTAGAAGAGTTAGGTATCAATTTTGA  
CCACATCTGGCAAAAAACACTAACAGTATTACATCCGTTAAAAGTAGCAGATGGCAGCATCATGAATGAA  
ACTGATTTGGCAGGTCCAATGGTTTTTGCCTTGCTTTTGGAGCCACATTGCTACTGGCTGGCAAAATCC  
AGTTTGGCTATGTATACGGGATCAGTGAATTGGATGTCTAGGAATGTTTTGTTTAACTTAATGAG  
TATGACAGGTGTTTCATTTGGTTGTGTGCAAGTGTCTTGGATATTGTCTTCTGCCCATGATCCTACTT  
TCCAGCTTTGAGTGATATTTCTTTGCAAGGAATGGTAGGAATCATTCTCACTGCTGGGATTATTGGAT  
GGTGTAGTTTTCTGCTTCCAAAATATTTCTGCATTAGCCATGGAAGGACAGCAACTTTTAGTAGC  
ATATCCTTGGCCTTTGTTATATGGAGTCTTGCCTGATTTCCGTCTTT

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC232110 representing NM\_001271732  
**Red=Cloning site Green=Tags(s)**

MMQPQQPYTGQIYQPTQAYTPASPQPFYGNFDEPPLLEELGINFDHIWQKTLTVLHPLKVADGSIMNE  
TDLAGPMVFLAFGATLLLAGKIQFGYVYGISAIIGCLGMFCLLNMSMTGVSFGCVASVLGYCLLPMILL  
SSFVAVIFSLQGMVGIILTAGIIGWCSFSASKIFISALAMEGQQLLVAYPCALLYGVFALISVF

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001271732

ORF Size: 609 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001271732.2](#)

**RefSeq Size:** 3343 bp

**RefSeq ORF:** 612 bp

**Locus ID:** 81555

**UniProt ID:** [Q969M3](#)

**Cytogenetics:** 5q31.3

**Protein Families:** Transmembrane

**MW:** 22.3 kDa

**Gene Summary:** Plays a role in transport between endoplasmic reticulum and Golgi.[UniProtKB/Swiss-Prot Function]