

Product datasheet for RC232126

AAGAB (NM_001271886) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: AAGAB (NM_001271886) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: AAGAB
Synonyms: KPPP1; p34; PPKP1; PPKP1A
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232126 representing NM_001271886
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATCTTGGTCTGCGATAGAGTGTCTGAAGATGGTATAAACCGACAAAAAGCTCAAGAATGGTGCATCA
AACATGGCTTTGAATTGGTAGAACTTAGTCCAGAGGAGTTGCCTGAGGAGGATGATGACTCCCAGAATC
TACAGGAGTAAAGCGAATTGTCCAAGCCCTGAATGCCAATGTGTGGTCCAATGTAGTGATGAAGAATGAT
AGGAACCAAGGCTTTAGCCTTCTCAACTCATTGACTGGAACAAACCATAGCATTGGGTGAGCAGATCCCT
GTCACCCAGAGCAACCCATTTGCCAGCAGCAGATAGTACTGAATCCCTCTCTGATCATCGGGTGGTGC
ATCTAACACAACAGATGCCAGGTTGATAGCATTGTGGATCCCATGTTAGATCTGGATATTCAAGAATTA
GCCAGTCTTACCACTGGAGGAGGAGATGTGGAGAAATTTGAAAGACTCTTTTCAAAGTTAAAGGAAATGA
AAGACAAGGCTGCGACGCTTCTCATGAGCAAAGAAAAGTGCATGCAGAAAAGGTGGCCAAAGCATTCTG
GATGGCAATCGGGGGAGACAGAGATGAAATTGAAGGCCTTTCATCTGATGAAGAGCAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC232126 representing NM_001271886
Red=Cloning site Green=Tags(s)

MILVCDRVSIEDGINRQKAQEWCIKHGFELVELSPEELPEEDDDFPESTGVKRIQALNANVSNVVMKND
RNQGFSLNLSLTGTNHSIGSADPCHPEQPHLPAADSTESLSDHRGGASNTTDAQVDSIVDPLDLDIQEL
ASLTTGGDVENFERLF SKLKEMKDKAATLPHEQRKVHAEKVAKAFWMAIGGDRDEIEGLSSDEEH


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:

Plasmid Map:


ACCN: NM_001271886

ORF Size: 618 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001271886.1](#), [NP_001258815.1](#)

RefSeq Size: 3344 bp

RefSeq ORF: 621 bp

Locus ID: 79719

UniProt ID: [Q6PD74](#)

Cytogenetics: 15q23

MW: 23.2 kDa

Gene Summary: The protein encoded by this gene interacts with the gamma-adaptin and alpha-adaptin subunits of complexes involved in clathrin-coated vesicle trafficking. Mutations in this gene are associated with type I punctate palmoplantar keratoderma. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Dec 2012]