

Product datasheet for RC232139

Geminin (GMNN) (NM_001251991) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Geminin (GMNN) (NM_001251991) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Geminin
Synonyms:	Gem; MGORS6
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC232139 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATCCAGTATGAAGCAGAAACAAGAAGAAATCAAAGAGAATATAAAGAATAGTTCTGTCCCAAGAA
GAACTCTGAAGATGATTCAGCCTTCTGCATCTGGATCTCTTGTGGAGAGAAAAAGAGCTGTCCGAGG
CTTGTCCAAAAGGAAACATCGGAATGACCACTAACATCTACAACCTCCAGCCCTGGGGTTATTGTCCCA
GAATCTAGTAAAATAAAAAATCTTGGAGGAGTCACCCAGGAGTCATTTGATCTTATGATTAAGAAAAATC
CATCCTCTCAGTATTGGAAGGAAGTGGCAGAAAAACGGAGAAAGGCGCTGTATGAAGCACTTAAGGAAAA
TGAGAAATTCATAAAGAAATTGAACAAAAGGACAATGAAATTGCCCGCCTGAAAAAGGAGAATAAAGAA
CTGGCAGAAGTAGCAGAACATGTACAGTATATGGCAGAGCTAATAGAGAGACTGAATGGTGAACCTCTGG
ATAATTTTGAATCACTGGATAATCAGGAATTTGATTCTGAAGAAGAAACTGTTGAGGATTCTAGTGGAA
AGACTCAGAAATTGGCACGTGTGCTGAAGGAAGTATCTTCTCTACGGATGCAAAGCCATGTATA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>RC232139 protein sequence Red=Cloning site Green=Tags(s)
-------------------	--

MNPSMKQKQEEIKENIKNSSVPRRTLKMIQPSASGSLVGRENELSAGLSKRKHRNDHLTSTTSSPGVIVP
ESSENKNLGGVTQESFDLMIKENPSSQYWKEVAEKRRKALYEALKENEKLHKEIEQKDNEIARLKKENKE
LAEVAEHVQYMAELIERLNGEPLDNFESLDNQEFDSEETTVEDSLVEDSEIGTCAEGTVSSSTDAKPCI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Chromatograms: https://cdn.origene.com/chromatograms/mk6145_b05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001251991

ORF Size: 627 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001251991.1](#), [NP_001238920.1](#)

RefSeq Size: 1051 bp


RefSeq ORF: 630 bp

Locus ID: 51053

UniProt ID: [O75496](#)
Cytogenetics: 6p22.3
Protein Families: Druggable Genome, Stem cell - Pluripotency
MW: 23.6 kDa

Gene Summary: This gene encodes a protein that plays a critical role in cell cycle regulation. The encoded protein inhibits DNA replication by binding to DNA replication factor Cdt1, preventing the incorporation of minichromosome maintenance proteins into the pre-replication complex. The encoded protein is expressed during the S and G2 phases of the cell cycle and is degraded by the anaphase-promoting complex during the metaphase-anaphase transition. Increased expression of this gene may play a role in several malignancies including colon, rectal and breast cancer. Alternatively spliced transcript variants have been observed for this gene, and two pseudogenes of this gene are located on the short arm of chromosome 16. [provided by RefSeq, Oct 2011]

Product images:



Circular map for RC232139