

Product datasheet for **RC232156**

ZFAND2B (NM_001270999) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ZFAND2B (NM_001270999) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ZFAND2B
Synonyms: AIRAPL
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232156 representing NM_001270999
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTTCCGGACCTCGGCGCTCACTGTTCCGAGCCGAGCTGTCAGCGCTTGGATTTCTGCCGCTTA
AGTGTGATGCCTGCTCAGGCATCTTCTGCGCAGACCATGTGGCCTACGCCCAGCATCACTGTGGATCTGC
TTACAAAAGATCTTACCAATAAGTGTGAACGCGCTGGCTGCCGCGAGGAGAAATGATGAACTGACC
TGTGAACGCTGTAGCCGAACTTCTGCATCAAGCACCGGCATCCACTGGACCATGATTGCTCTGGGGAGG
GGCACCCAACAGCCGGGAGGACTTGCTGCCATCTCCAGAGCACAAGCTGTGGCTTCTACAAGCACTGT
CCCCAGCCCAAGTCAAACCATGCCTTCTGTACCTCTCCAGCAGAGCCACAACCCGATCTCCGTCTCGG
ACAGCCCCCTCAGTGATTGCTTTGCAGAATGGCCTGAGTGAGGATGAAGCTCTGCAGCGGGCCCTGGAAA
TGTCCCTGGCAGAAACCAACCCAGGTTCCAAGTTGTCAGGAGGAAGAAGACCTAGCTTTAGCACAAAGC
ACTGTCAGCCAGTGAGGCAGAATACCAGCGGCAGCAGGCCAGAGCCGAGCTCGAAGCCGTCCAAGTGC
AGCCTGTGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC232156 representing NM_001270999
 Red=Cloning site Green=Tags(s)

MEFPDLGAHCSEPCQRLDFLPLKCDACSGIFCADHVAYAQHHCGSAYQKIFTNKCRAGCRQREMMKLT
 CERCSRNF CIKRRHPLDHDCSGEGHPTSRAGLA AISRAQAVASTVPSPTSQTMPSTSPSRATTRPSW
 TAPPVIALQNLSEDEALQRALEMSLAETKPQVPSCQEEEDLALAQALSASEAEYQRQAQSRSSKPSNC
 SLC


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001270999

ORF Size: 639 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001270999.1 , NP_001257928.1
RefSeq Size:	1070 bp
RefSeq ORF:	642 bp
Locus ID:	130617
UniProt ID:	Q8WV99
Cytogenetics:	2q35
MW:	23.6 kDa
Gene Summary:	This gene encodes a protein containing AN1-type zinc-fingers and ubiquitin-interacting motifs. The encoded protein likely associates with the proteasome to stimulate the degradation of toxic or misfolded proteins. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Aug 2012]