

Product datasheet for **RC232182**

ARD1A (NAA10) (NM_001256119) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ARD1A (NAA10) (NM_001256119) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: NAA10
Synonyms: ARD1; ARD1A; ARD1P; DXS707; hARD1; MCOPS1; NATD; OGDNS; TE2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232182 representing NM_001256119
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAACATCCGCAATGCGAGGCCAGAGGACCTAATGAACATGCAGCACTGCAACCTCCTCTGCCTGCCCG
AGAACTACCAGATGAAATACTACTTCTACCATGGCCTTCTCTGGCCCCAGCTCTTTACATTGCTGAGGA
CGAGAATGGGAAGATTGTGGGGTATGTCCTGGCCAAAATGGAAGAGGACCCAGATGATGTGCCCATGGA
CATATCACCTCATTGGCTGTGAAGCGTTCCACCGGCCTCGGTCTGGCTCAGAACTGATGGACCAGG
CCTCTCGAGCCATGATAGAGAATTCAATGCCAAATATGTCTCCCTGCATGTCAGGAAGAGGATCAGTGA
AGTGGAGCCAAATACTATGCAGATGGGGAGGACGCCTATGCCATGAAGCGGGACCTCACTCAGATGGCC
GACGAGCTGAGGCGGCACCTGGAGCTGAAAGAGAAGGGCAGGCACGTGGTGGTGGTGGCCATCGAGAACA
AGGTGGAGAGCAAAGCAATTCACCTCCGAGCTCAGGAGAGGCCGTGTCGCGAGGAGAAGGGCCTGGCTGC
CGAGGATAGTGGTGGGGACAGCAAGGACCTCAGCGAGGTCAGCGAGACCACAGAGACACAGATGTCAAG
GACAGCTCAGAGGCTCCGACTCAGCCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC232182 representing NM_001256119
 Red=Cloning site Green=Tags(s)

MNIRNARPEDLMNQHCNLLCLPENYQMKYYFYHGLSWPQLSYIAEDENKIVGVVLAKMEEDPDDVPHG
 HITSLAVKRSHRRLGLAQKLMQASRAMIENFNAKYVSLHVRKRISEVPEPKYYADGEDAYAMKRDLTQMA
 DELRRHLELKEKGRHVVLGAIENKVESKGNPPSSGGEACREEKGLAAEDSGGDSKDLSEVSETTESTDVK
 DSSEASDSAS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001256119

ORF Size: 660 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256119.2
RefSeq Size:	1091 bp
RefSeq ORF:	663 bp
Locus ID:	8260
UniProt ID:	P41227
Cytogenetics:	Xq28
Protein Families:	Druggable Genome
Protein Pathways:	Glycerophospholipid metabolism, Limonene and pinene degradation, Phenylalanine metabolism, Tyrosine metabolism
MW:	25.2 kDa
Gene Summary:	N-alpha-acetylation is among the most common post-translational protein modifications in eukaryotic cells. This process involves the transfer of an acetyl group from acetyl-coenzyme A to the alpha-amino group on a nascent polypeptide and is essential for normal cell function. This gene encodes an N-terminal acetyltransferase that functions as the catalytic subunit of the major amino-terminal acetyltransferase A complex. Mutations in this gene are the cause of Ogden syndrome. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Jan 2012]