

## Product datasheet for **RC232190**

### **TMEM235 (NM\_001204210) Human Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** TMEM235 (NM\_001204210) Human Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** TMEM235  
**Synonyms:** ARGM1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RC232190 representing NM\_001204210  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCCCGCTGGGCGCGCTGCTCCTGGCCGCCCCCTGGGTGCACTGCTCAGCTTCGGCTCCTGGCCG  
CCGCGGTCGCCAGCGACTACTGGTACATCCTGGAGGTGGCGGACCCGGCAATGGCAGCGCTGGCCCG  
GCGCGCAGAGCTGCTCCTCGCACTCGGGGCTCTGGCGCATCTGCGAAGGGCAGAACGGCTGCATCCCG  
CTGGTCGACCCTTTTCCAGTGAGAGCCTGGACGTCTCCACCTCGGTGCAGCACCTCATCTTGCTGCACC  
GTGCAGTCATTGTGGTCCTGCCCTGAGCCTGGTCTTCTCGTGTGTGGCTGGATCTGCGCCTGCTCAG  
CTCCCTGGCCAGAGCGTGTCTGCTGCTTTTACCAGGCTGCTACTTCTGCTGGGAGTGTCTGACA  
CTGGCGGGGGTACAGATCTACATCAGCTACTCGCACCTGGCCTTTGCGGAGACGGTGCAGCAGTATGGCC  
CGCAGCACATGCAGGGCGTCCGCGTCAGCTTCCGCTGGTCCATGGCCCTGGCCTGGGGCTCCTGTGCCTT  
GGAGGCATTACAGCGGAACCCTCCTGCTCTCAGCTGCCTGGACCCCTCAGCCTGAGCCCCCAATCTGTGGT  
CATCTGAGTCCCCAGCAGGTGGGAGGGAGAGGGGAGAC

**ACGCGT**ACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC232190 representing NM\_001204210  
Red=Cloning site Green=Tags(s)

MARLGALLLAAALGALLSFALLAAAVASDYWYILEVADAGNGSAWPGRAELLSHSGLWRICEGQNGCIP  
 LVDPFAESLDVSTSVQHLILLHRAVIVLPLSLVLLVCGWICGLLSSLAQSVLLLLFTGCYFLLGSVLT  
 LAGVSIYISYSHLAFETVQYGPQHMQGVVRSFGWSMALAWGSCALEAFSGTLLLSAAWTLSSLSPPICG  
 HLLSPQQVGGRRGGD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001204210

**ORF Size:** 669 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001204210.1</a> , <a href="#">NP_001191139.1</a>
<b>RefSeq Size:</b>	2509 bp
<b>RefSeq ORF:</b>	672 bp
<b>Locus ID:</b>	283999
<b>UniProt ID:</b>	<a href="#">A6NFC5</a>
<b>Cytogenetics:</b>	17q25.3
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	23.8 kDa