

Product datasheet for **RC232208**

PAMM (FAM213A) (NM_001243778) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PAMM (FAM213A) (NM_001243778) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PRXL2A
Synonyms: Adrx; C10orf58; FAM213A; PAMM
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232208 representing NM_001243778
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCTTCTCCAGGACCAAGTTTCTCACCATGGGGATGTGGTCCATTGGTGCAGGAGCCCTGGGG
CTGCTGCCTTGGCATTGCTGCTTGCCAACACAGACGTGTTTCTGTCCAAGCCCCAGAAAGCGGCCCTGGA
GTACCTGGAGGATATAGACCTGAAAACACTGGAGAAGGAACCAAGGACTTCAAAGCAAAGGAGCTATGG
GAAAAAATGGAGCTGTGATTATGGCCGTGCGGAGGCCAGGCTGTTTCTCTGTGAGAGGAAGCTGCGG
ATCTGTCTCCCTGAAAAGCATGTTGGACCAGCTGGGCGTCCCCCTCTATGCAGTGGTAAAGGAGCACAT
CAGGACTGAAGTGAAGGATTTCCAGCCTTATTTCAAAGGAGAAATCTTCCTGGATGAAAAGAAAAGTTC
TATGGTCCACAAAGGCGGAAGATGATGTTTATGGGATTTATCCGTCTGGGAGTGTGGTACAACCTTCTCC
GAGCCTGGAACGGAGGCTTCTCTGGAACCTGGAAGGAGAAGGCTTCATCCTTGGGGGAGTTTTCTGGT
GGGATCAGGAAAGCAGGGCATTCTTCTGAGCACCGAGAAAAAGAATTTGGAGACAAAGTAAACCTACTT
TCTGTTCTGGAAGCTGCTAAGATGATCAAACCACAGACTTGGCCTCAGAGAAAAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001243778.1 , NP_001230707.1
RefSeq Size:	1673 bp
RefSeq ORF:	690 bp
Locus ID:	84293
UniProt ID:	Q9BRX8
Cytogenetics:	10q23.1
Protein Families:	Secreted Protein, Transmembrane
MW:	26.2 kDa
Gene Summary:	Involved in redox regulation of the cell (PubMed:26438880, PubMed:19951071). Acts as an antioxidant (PubMed:19951071, PubMed:26438880). Inhibits TNFSF11-induced NFkB1 and JUN activation and osteoclast differentiation (PubMed:19951071). May affect bone resorption and help to maintain bone mass (PubMed:19951071). Acts as a negative regulator of macrophage-mediated inflammation by inhibiting macrophage production of inflammatory cytokines, probably through suppression of the MAPK signaling pathway (PubMed:26438880). [UniProtKB/Swiss-Prot Function]