

Product datasheet for RC232221

CHODL (NM_001204176) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CHODL (NM_001204176) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	CHODL
Synonyms:	C21orf68; MT75; PRED12
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC232221 representing NM_001204176 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCTACTTCCATGAAGTGTCCAGCCGAGTGAGCTTTCAGGAGGCACGCCTGGCTTGTGAGAGTGAGG
GAGGAGTCCTCCTCAGCCTTGAGAATGAAGCAGAACAGAAGTTAATAGAGAGCATGTTGCAAACCTGAC
AAAACCCGGGACAGGGATTTCTGATGGTGATTTCTGGATAGGGCTTTGGAGGAATGGAGATGGGCAAACA
TCTGGTGCCTGCCAGATCTCTACCAGTGGTCTGATGGAAGCAATCCCAAGTACCGAAACTGGTACACAG
ATGAACCTTCTGCGGAAGTAAAAGTGTGTTGTGATGTATCACCAACCAACTGCCAATCCTGGCCTTGG
GGGTCCCTACCTTACCAGTGAATGATGACAGGTGTAACATGAAGCACAATTATTTGCAAGTATGAA
CCAGAGATTAATCCAACAGCCCTGTAGAAAAGCCTTATCTTACAAATCAACCAGGAGACACCCATCAGA
ATGTGGTTGTTACTGAAGCAGGTATAATTCCCAATCTAATTTATGTTGTTATACCAACAATACCCCTGCT
CTTACTGATACTGGTTGCTTTTGGAACTGTTGTTCCAGATGCTGCATAAAAAGTAAAGGAAGAACAAAA
ACTAGTCCAAACCAGTCTACACTGTGGATTTCAAAGAGTACCAGAAAAGAAAGTGGCATGGAAGTA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC232221 representing NM_001204176
 Red=Cloning site Green=Tags(s)

MAYFHELSSRVSFQEARLACESEGGVLLSLENAEQKLIESMLQNLTKPGTGISDGFWIGLWRNGDGQT
 SGACPDLYQWSDGSNSQYRNWYTDEPSCGSEKCVVMYHQPTANPGLGGPYLYQWDDRCNMKHNYICKYE
 PEINPTAPVEKPYLTNQPGDTHQNVVVTEAGIIPNLIYVVIPTIPLLLLILVAFGTCCFQMLHKSGRTK
 TSPNQSTLWISKSTRKESGMEV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_001204176

ORF Size: 696 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001204176.1 , NP_001191105.1
RefSeq Size:	2316 bp
RefSeq ORF:	699 bp
Locus ID:	140578
UniProt ID:	Q9H9P2
Cytogenetics:	21q21.1
Protein Families:	Transmembrane
MW:	26.4 kDa
Gene Summary:	This gene encodes a type I membrane protein with a carbohydrate recognition domain characteristic of C-type lectins in its extracellular portion. In other proteins, this domain is involved in endocytosis of glycoproteins and exogenous sugar-bearing pathogens. This protein localizes predominantly to the perinuclear region. Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Feb 2011]