

Product datasheet for **RC232230**

RAB18 (NM_001256410) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RAB18 (NM_001256410) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: RAB18
Synonyms: RAB18LI1; WARBM3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232230 representing NM_001256410
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACGAGGACGTGCTAACCACCCTGAAGATCCTCATCATCGGCGAGAGTGGGGTGGCAAGTCCAGCC
TGCTCTTGAGGTTACAGATGATACGTTTGATCCAGAACTGCAGCAACAATAGGTGTTGACTTTAAGGT
GAAAACAATTTCACTGGATGAAATAAGGCTAACTTGAATATGGGTTACTTTGCATCAGCAGACTGCA
AACTTTTTCTAAAGAGCCAGATAGGAAATAGTCCAATCCTGAAGTGGCCATGTGGCAATACGATACTG
CTGGTCAAGAGAGGTTTAGAACATTAACCTCCAGCTATTATAGAGGTGCACAGGGTGTATATTAGTTTA
TGATGTACAAGAAGAGATACATTTGTTAACTGGATAATTGGTTAAATGAATTGGAAACATACTGTACA
AGAAATGACATAGTAAACATGCTAGTTGGAAATAAAATCGATAAGGAAAAATCGTGAAGTCGATAGAAATG
AAGGCCTGAAATTTGCACGAAAGCATTCCATGTTATTTATAGAGGCAAGTGCAAAAACCTGTGATGGTGT
ACAATGTGCCTTTGAAGAACTTGTGAAAAGATCATTAGACCCCTGGACTGTGGGAAAGTGAGAACCAG
AATAAAGGAGTCAAATGTCACACAGGGAAGAAGGCCAAGGAGGAGCCTGTGGTGGTTATTGCTCTG
TGTTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RC232230 representing NM_001256410
 Red=Cloning site Green=Tags(s)

MDEDVLTTLKILIIIGESGVGKSSLLLRFTDDTFDEPELAATIGVDFKVKTISVDGNKAKLAIWVTLHQQT
 NFFLKSQIGNSPIKAWMWQYDTAGQERFRTLTPSYRGAQGVILVYDVTRRDTFVKLDNWLNELEYCT
 RNDIVNMLVGNKIDKENREVDREGLKFARKHSMLEASAKTCDGVQCAFEELVEKIIQTPGLWESENQ
 NKGVKLSHREEGQGGACGGYCSVL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001256410
ORF Size: 705 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001256410.2
RefSeq Size:	5090 bp
RefSeq ORF:	708 bp
Locus ID:	22931
UniProt ID:	Q9NP72
Cytogenetics:	10p12.1
Protein Families:	Druggable Genome
MW:	26.9 kDa
Gene Summary:	The protein encoded by this gene is a member of a family of Ras-related small GTPases that regulate membrane trafficking in organelles and transport vesicles. Knockdown studies in zebrafish suggest that this protein may have a role in eye and brain development. Mutations in this gene are associated with Warburg micro syndrome type 3. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jan 2012]