

Product datasheet for **RC232406**

GSG1 (NM_001206843) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: GSG1 (NM_001206843) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: GSG1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232406 representing NM_001206843
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCCAAGATGGAGCTCTCGAAGGCCCTTCTTGCCAGCGGACACTCCTATCTGCCATCCTCAGCATGC
TATCACTCAGCTTCTCCACAACATCCCTGCTCAGCAACTACTGGTTTGTGGGCACACAGAAGGTGCCCAA
GCCCTGTGCGAGAAAGGTCTGGCAGCCAAGTGCTTTGACATGCCAGTGTCCCTGGATGGAGATACCAAC
ACATCCACCCAGGAGGTGGTACAATACTGGGAGACTGGGGATGACCGGTTCTCCTCCGGAGCTTCC
GGAGTGGCATGTGGCTATCCTGTGAGGAACTGTGGAAGAACCAGCACTGCTCCATCCCCAGTCCTGGAA
ACAATTTAGAGCCCTTCGGTCCAGTGGTACAGCGGCAGCAAAAGGGGAGAGGTGCCGAAGTTTCATTGAA
CTTACACCACCAGCCAAGAGAGGTCTCCTGGGGATGGTGGCCACATGATGTATTACAAGTCTTCCAAG
CGACTGTCAACTTGGGTCCAGAAGACTGGAGACCACATGTTTGAATTATGGCTGGGCTTCTACATGGC
CTGGCTCCTTACCTGCTGCATGGCGTCGGCTGTCAACACGTACACCAGGATGGTGTG
GAGTTCAAGTCAAGCATAGTAAGAGCTTCAAGGAAAACCGAACTGCCTACCACATCACCATCAGTGTT
TCCCTCGGCGGCTGTCAAGTGCAGCCCCACCGTGGGTCTTTGACCAGTACCACCAGTATCATAATCA
GCCCATCCACTCTGTCTGAGGGAGTCGACTTCTACTCCGAGCTGCGGAACAAGGGATTTCAAAGAGGG
GCCAGCCAGGAGCTGAAAGAAGCAGTTAGGTCATCTGTAGAGGAAGAGCAGTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

ORF Size:	894 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001206843.3
RefSeq Size:	2460 bp
RefSeq ORF:	897 bp
Locus ID:	83445
UniProt ID:	Q2KHT4
Cytogenetics:	12p13.1
Protein Families:	Transmembrane
MW:	34.3 kDa
Gene Summary:	May cause the redistribution of PAPOLB from the cytosol to the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function]