

Product datasheet for **RC232431**

ARFIP2 (NM_001242855) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ARFIP2 (NM_001242855) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ARFIP2
Synonyms: POR1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232431 representing NM_001242855
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTGTCAAGGACCACTCAATGAAACCAGCATTGTGTCTGGTGGCTATGGGGCTCTGGTGATGGAC
TCATCCCCACAGGGTCTGGCCGCATCCATCTCACAGCACCCTCTTCTGGCCCTGGAGATGAGGTGGC
TCGGGGCATTGCTGGAGAAAAGTTTGACATCGTCAAGAAATGGGGCATCAACACCTATAAGTGCACAAAG
CAACTGTTATCAGAACGATTTGGTCGAGGCTCACGGACTGTGGACCTGGAGCTAGAGCTGCAGATTGAGT
TGCTGCGTGAGACGAAGCGCAAGTATGAGAGTGTCTGCAGCTGGGCCGGGCACTGACAGCCACCTCTA
CAGCCTGCTGCAGACCCAGCATGCACTGGGTGATGCCTTTGCTGACCTCAGCCAGAAGTCCCCAGAGCTT
CAGGAGGAATTTGGCTACAATGCAGAGACACAGAACTACTATGCAAGAATGGGGAAACGCTGCTAGGAG
CCGTGAACCTCTTTGCTCTAGCATCAACACATTGGTCACCAAGACCATGGAAGACACGCTCATGACTGT
GAAACAGTATGAGGCTGCCAGGCTGGAATATGATGCCTACCGAACAGACTTAGAGGAGCTGAGTCTAGGC
CCCCGGGATGCAGGGACACGTGGTCGACTTGAGAGTGCCAGGCCACTTCCAGGCCCATCGGGACAAGT
ATGAGAAGCTGCGGGGAGATGTGGCCATCAAGCTCAAGTTCTTGGAAAGAAAACAAGATCAAGGTGATGCA
CAAGCAGCTGCTGCTTCCACAATGCTGTGTCCGCTACTTTGCTGGGAACCAGAAACAGCTGGAGCAG
ACCTGCAGCAGTTCAACATCAAGCTGCGGCCTCCAGGAGCTGAGAAACCTCCTGGCTAGAGGAGCAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC232431 representing NM_001242855
 Red=Cloning site Green=Tags(s)


MVSGPNLNETSIVSGGYGGSGDGLIPTGSGRHPSSHSTTPSGPGDEVARGIAGEKFDIVKKWGINTYKCTK
 QLLSERFGRGSRVTDLELELQIELLRETKRKYESVLQLGRALTAHLYSLLQTQHALGDAFADLSQKSP
 QEEFGYNAETQKLLCKNGETLLGAVNFFVSSINTLVTKTMDTLMTVKQYEAARLEYDAYRTDLEELSLG
 PRDAGTRGRLESAQATFQHRDKYEKLRGDVAIKLKFLEENKIKVMHKQLLLL FHNVAVSAYFAGNQKQLEQ
 TLQQFNIKL RPPGAEKPSWLEEQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001242855

ORF Size: 909 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001242855.3
RefSeq Size:	3624 bp
RefSeq ORF:	912 bp
Locus ID:	23647
UniProt ID:	P53365
Cytogenetics:	11p15.4
Protein Families:	Druggable Genome
MW:	34.2 kDa
Gene Summary:	Putative target protein of ADP-ribosylation factor. Involved in membrane ruffling. [UniProtKB/Swiss-Prot Function]