

Product datasheet for **RC232440**

ADK (NM_001202450) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ADK (NM_001202450) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ADK
Synonyms: AK
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232440 representing NM_001202450
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGCTGCTGAGGAGGAGCCGAAGCCAAAAAGCTGAAGGTGGAGGCGCCGAAGCGCTGAGAGAAA
ATATTCTCTTTGGAATGGGAAATCCTCTGCTTGACATCTCTGCTGTAGTGGACAAAGATTCCTTGATAA
GTATTCTCTGAAACCAAATGACCAAATCTGGCTGAAGACAAACACAAGGAAGTGTGATGAAGTGTG
AAAAAATTCAAAGTCGAATATCATGCTGGTGGCTCTACCCAGAATTCAAATTAAGTGGCTCAGTGGATGA
TTCAACAGCCACACAAAGCAGCAACATTTTTGGATGCATTGGGATAGATAAATTTGGGGAGATCCTGAA
GAGAAAAGCTGCTGAAGCCCATGTGGATGCTCATTACTACGAGCAGAATGAGCAGCCAACAGGAAGTGT
GCTGCATGCATCACTGGTGACAACAGGTCCCTCATAGCTAATCTTGCTGCTGCCAATTGTTATAAAAAGG
AAAAACATCTTGATCTGGAGAAAACTGGATGTTGGTAGAAAAAGCAAGAGTTTGTATATAGCAGAAGC
TGCCACTTTTGTAGAGAGCAAGGCTTTGAGACTAAAGACATTAAGAGATAGCCAAAAAGACACAAGCC
CTGCCAAAGATGAACTCAAAGAGGCAGCGAATCGTGATCTCACCAAGGGAGAGATGACACTATAATGG
CTACAGAAAGTGAAGTCACTGCTTTTGGTGTCTGGATCAAGACCAGAAAGAAATTATTGATACCAATGG
AGCTGGAGATGCATTTGTTGGAGTTTTCTGTCTCAACTGGTCTCTGACAAGCCTCTGACTGAATGTATC
CGTGCTGGCCACTATGCAGCAAGCATCATAATTAGACGGACTGGCTGCACCTTCTGAGAAGCCAGACT
TCCAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC232440 representing NM_001202450
 Red=Cloning site Green=Tags(s)


MAAAEEEPKPKKLKVEAPQALRENILFGMGNPLLDISAVVDKDFLDKYSCLKPNDQILAEDKHKELFDELV
 KFKFVEYHAGGSTQNSIKVAQWMIQQPHKAATFFGCIGIDKFGEILKRKAAEAHVDAHYEQNEQPTGTC
 AACITGDNRSLIANLAAANCYKKEKHLDEKNWMLVEKARVCYIAEAATFAREQGFETKDIKEIAKKTQA
 LPKMNSKRQRIVIFTQGRDDTIMATESEVTAFAVLDDQDQKEIIDTNGAGDAFVGGFLSQLVSDKPLTECI
 RAGHYAASIIIRRTGCTFPEKPDFH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001202450

ORF Size: 915 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001202450.2
RefSeq Size:	1878 bp
RefSeq ORF:	918 bp
Locus ID:	132
UniProt ID:	P55263
Cytogenetics:	10q22.2 10q11-q24
Protein Families:	Druggable Genome
Protein Pathways:	Metabolic pathways, Purine metabolism
MW:	34.5 kDa
Gene Summary:	This gene an enzyme which catalyzes the transfer of the gamma-phosphate from ATP to adenosine, thereby serving as a regulator of concentrations of both extracellular adenosine and intracellular adenine nucleotides. Adenosine has widespread effects on the cardiovascular, nervous, respiratory, and immune systems and inhibitors of the enzyme could play an important pharmacological role in increasing intravascular adenosine concentrations and acting as anti-inflammatory agents. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2011]