

Product datasheet for RC232479

Natriuretic Peptide Receptor C (NPR3) (NM_001204376) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Myc-DDK
Symbol:	Natriuretic Peptide Receptor C
Synonyms:	ANP-C; ANPR-C; ANPRC; C5orf23; GUCY2B; NPR-C; NPRC
Mammalian Cell	Neomycin
Selection:	
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)

ORF Nucleotide Sequence: >RC232479 representing NM_001204376
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCTGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGCCATCTGCACTGGGACCTTGGTCCTTGTTCCTGACCTTGCGCCGAGGGACCAGGAAGGTCAGG
TGCTCGCCCCGCGTCGGTGCTTCAGGAGGGAATGTGCCAGAACGAGAAAGTGGTGATCATGTGTGCGAG
CAGTGACACCATCCGGAGCATCATGCTGGTGCGCACAGGCATGGCATGACCACTGGAGACTACGCCTTC
TTCAACATTGAGCTCTTCAACAGCTCTTCTATGGAGATGGCTCATGGAAGAGAGGAGACAAACAGCACT
TTGAAGCTAAGCAAGCACTACTCGTCCCTCCAGACAGTCACTCTACTGAGGACAGTGAAACCTGAGTTTGA
GAAGTTTTCCATGGAGGTGAAAAGTTCAAGTTGAGAAAACAAGGGCTCAATATGGAGGATTACGTTAATG
TTTGTGAAGGATTCCACGATGCCATCCTCCTCTACGTCTTGGCTCTACATGAAGTACTCAGAGCTGGTT
ACAGCAAAAAGGATGGAGGGAAAATTATACAGCAGACTTGGAACAGAACATTTGAAGGTATCGCCGGGCA
GGTGTCCATAGATGCCAACGGAGACCGATATGGGGATTTCTCTGTGATTGCCATGACTGATGTGGAGGCG
GGCACCCAGGAGGTTATTGGTGATTATTTTGGAAAAGAAGGTCGTTTTGAAATGCGGCCGAATGTCAAAT
ATCCTTGGGGCCCTTTAAACTGAGAATAGATGAAAACCGAATTGTAGAGCATACAAACAGCTCTCCCTG
CAAATCATGTGGCCTAGAAGATCGGCAGTGACAGGAATTGTCGTGGGGCTTTACTAGGAGCTGGCTTG
CTAATGGCCTTCTACTTTTTCAGGAAGAAATACAGAAATACCAATTGAGAGGCGAACCCAGCAAGAAGAAA
GTAACCTTGAAAACATCGGGAATTACGGGAAGATTCCATCAGATCCCATTTCAGTAGCT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence:

>RC232479 representing NM_001204376

Red=Cloning site Green=Tags(s)


MEPSALGPWSLFPDLAPRDQEGQVLAPRRCFRREGQNEKVVIMCASSDTIRSIMLVHRHGMTSGDYAF
FNIELFNSSSYGDGSKRGDKHDFEAKQAYSSLQTVTLRLTVKPEFEKFSMEVKSSVEKQGLNMEDYVNM
FVEGFHDAILLVYALHEVL RAGYSKDGKIIQQTWNRTFEGIAQVSDANGDRYGD FSVIAMTDVEA
GTQEVIGDYFGKEGRFEMRPNVKYPWGPLKL RIDENRIVEHTNSSPCKSCGLEESAVTGIVVGALLGAGL
LMAFYFFRKKYRITIERRTQQEESNLGKHRELREDSIRSHFSVA

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN:

NM_001204376

ORF Size:

972 bp

OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:

This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:

The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001204376.1](#), [NP_001191305.1](#)

RefSeq Size: 6380 bp

RefSeq ORF: 975 bp

Locus ID: 4883

UniProt ID: [P17342](#)


Cytogenetics: 5p13.3

Protein Families: Druggable Genome, Transmembrane

MW: 37.1 kDa

Gene Summary: This gene encodes one of three natriuretic peptide receptors. Natriuretic peptides are small peptides which regulate blood volume and pressure, pulmonary hypertension, and cardiac function as well as some metabolic and growth processes. The product of this gene encodes a natriuretic peptide receptor responsible for clearing circulating and extracellular natriuretic peptides through endocytosis of the receptor. Multiple transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Feb 2011]

Product images:



Circular map for RC232479