

Product datasheet for **RC232493**

PRPSAP2 (NM_001243936) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: PRPSAP2 (NM_001243936) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: PRPSAP2
Synonyms: PAP41
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232493 representing NM_001243936
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTTTTGTGTGACGCCACCTGAATTAGAAACCAAGATGAACATAACCAAAGGTGGTCTGGTGTGTTTT
CAGCAAATCGAATTCATCATGTATGGAGCTATCAAAGAAAATTGCAGAGGACGTGAACACCACCATCAT
GGAGCTCCTGATCATGGTGTATGCATGTAAGACCTCTTGTGCCAAGAGCATCATTGGCGTGATACCCTCA
TTTCCTTACAGCAAGCAGTGCAAGATGAGAAAAAGAGGCTCCATTGTCTCTAAATTGCTGGCTTCCATGA
TGTGCAAAGCTGGTCTAACTCATCTTATTACTATGGATTTACACCAGAAGGAAATTCAGGGCTTCTTCAA
TATTCTGTGACAATTTAAGAGCATCTCCCTTCTTATTACAGTATATTCAAGAAGAGATCCCAGATTAC
AGGAATGCAGTAATCGTGGCCAAGTCTCCAGCCTCGGCGAAGAGGGCACAGTCTTTTGTGAGCGCCTGC
GCCTGGGAATTGCAGTGATTCATGGAGAGGCGCAGGATGCCGAGTCGGACTTGGTGGATGGACGGCATT
CCCACCCATGGTCAGAAGTGTGGCTGCCATCCACCCAGCCTGGAGATCCCCATGCTGATTCCTAAAGAA
AAGCCCCAATCACGGTGTGGGTGATGTTGGAGGAAGGATTGCCATCATCGTGGATGACATCATTGATG
ATGTTGACAGCTTCTTGTCTGCAGCAGAGACCCTGAAGGAAAGAGTGCATATAAGATCTTTGTGATGGC
AACTCATGGCTTGTGCTTCTGACGCCCCCGCGGATTGAAGAGTCTGCCATTGATGAGGTGGTGGTC
ACCAATACAATTCACATGAAGTCCAGAAGCTCCAGTGCCCAAGATTAAGTGTGGATATCAGCATGA
TCCTTTCAGAGGCGATCCGTCGGATCCACAATGGGGAGTCCATGTCTACCTTTTCAGAAACATAGGCTT
AGATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RC232493 representing NM_001243936
 Red=Cloning site Green=Tags(s)

MFCVTPPELETKMNITKGLVLFSA NSNSSCMELSKKIAEDVNTTIMELLIMVYACKTSCAKSIIGVIPY
 FPYSKQCKMRKRSIVSKLLASMMCKAGLTHLITMDLHQKEIQGFFNIPVDNLRASPFLQYIQEEIPDY
 RNAVIVAKSPASAKRAQSF AERLRLGIAVIHGEAQDAESDLVDGRHSPPMVRSVAAIHPSLEIPMLIPKE
 KPPITVVGDVGGRIAIIVDDIIDDVDSFLAAAE TLKERGAYKIFVMATHGLLSSDAPRRIEESAIDEVVV
 TNTIPHEVQKLQCPKIKTVDISMILSEAIRRHNGESMSYLFRNIGLDD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001243936

ORF Size: 987 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001243936.2
RefSeq Size:	1838 bp
RefSeq ORF:	990 bp
Locus ID:	5636
UniProt ID:	O60256
Cytogenetics:	17p11.2
Protein Families:	Druggable Genome
MW:	36.7 kDa
Gene Summary:	This gene encodes a protein that associates with the enzyme phosphoribosylpyrophosphate synthetase (PRS). PRS catalyzes the formation of phosphoribosylpyrophosphate which is a substrate for synthesis of purine and pyrimidine nucleotides, histidine, tryptophan and NAD. PRS exists as a complex with two catalytic subunits and two associated subunits. This gene encodes a non-catalytic associated subunit of PRS. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2011]