

Product datasheet for **RC232505**

ST3GAL4 (NM_001254758) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ST3GAL4 (NM_001254758) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ST3GAL4
Synonyms:	CGS23; gal-NAc6S; NANTA3; SAT3; SIAT4; SIAT4C; ST-4; ST3GalA.2; ST3GalIV; STZ
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC232505 representing NM_001254758 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGTCAGCAAGTCCCGCTGGAAGCTCCTGGCCATGTTGGCTCTGGTCTGGTCGTCATGGTGTGGTATT
CCATCTCCCGGAAGACAGGTACATCGAGCTTTTTTATTTCCCATCCCAGAGAAGAAGGAGCCGTGCC
CCAGGGTGAGGCAGAGAGCAAGGCCTCTAAGCTCTTTGGCAACTACTCCCGGGATCAGCCATCTTCCTG
CGGCTTGAGGATTATTTCTGGGTCAAGACGCCATCTGCTTACGAGCTGCCCTATGGGACCAAGGGGAGTG
AGGATCTGCTCCTCCGGGTGCTAGCCATCACCAGCTCCTCCATCCCAAGAACATCCAGAGCCTCAGGTG
CCGCCGCTGTGTGGTCTGTTGGGAACGGGACCGGCTGCGGAACAGCTCACTGGGAGATGCCATCAACAAG
TACGATGTGGTCATCAGATTGAACAATGCCCCAGTGGCTGGCTATGAGGGTGACGTGGGCTCCAAGACCA
CCATGCGTCTCTTCTACCCTGAATCTGCCCACTCGACCCAAAGTAGAAAAACAACCCAGACACTCCT
CGTCTGGTAGCTTTCAAGGCAATGGACTTCCACTGGATTGAGACCATCCTGAGTGATAAGAAGCGGGTG
CGAAAGGGTTCTGGAACAGCCTCCCTCATCTGGGATGTCAATCCTAAACAGATTTCGATTCTCAACC
CCTTCTCATGGAGATTGCAGCTGACAACTGCTGAGCCTGCAATGCAACAGCCACGGAAGATTAAAGCA
GAAGCCACACGGGCTGTTGGCCATCAGCTGGCCCTCCACCTCTGTGACTTGGTGACATTGCCGCGC
TTTGGCTACCCAGACGCCTACAACAAGAAGCAGACCATTCACTACTATGAGCAGATCAGCTCAAGTCCA
TGGCGGGGTGAGCCATAATGTCTCCCAAGAGGCCCTGGCCATTAAGCGGATGCTGGAGATGGGAGCTAT
CAAGAACCTCACGTCCTTC

ACGGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

Protein Sequence: >RC232505 representing NM_001254758
 Red=Cloning site Green=Tags(s)


MVSKSRWKL LAMLALVLVVMVWYSISREDRYIELFYFPIPEKKEPCLQGEAESKASKLFGNYSRDQPIFL
 RLEDYFVWKTPSAYELPYGKGSDDLRLVLAITSSSIPKNIQSLRCRRCVVVGNGHRLRNSSLGDAINK
 YDVVIRLNNAPVAGYEGDVGSKTTMRLFYPESAHFDPKVENNPDTLLVLVAFKAMDFHWIETILSDKKRV
 RKGFWKQPPLIWDVNPQIRILNPFMEIAADKLLSLPMQQPRKIKQKPTTGLLAI TLALHLCDLVHIAG
 FGYPDAYNKKQTIHYEQITLKS MAGSGHNVSQEALAIKRMLEMGAIKNLTSF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001254758

ORF Size: 999 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001254758.1 , NP_001241687.1
RefSeq Size:	1727 bp
RefSeq ORF:	1002 bp
Locus ID:	6484
UniProt ID:	Q11206
Cytogenetics:	11q24.2
Protein Families:	Secreted Protein, Transmembrane
Protein Pathways:	Glycosphingolipid biosynthesis - lacto and neolacto series, Metabolic pathways
MW:	38.5 kDa
Gene Summary:	This gene encodes a member of the glycosyltransferase 29 family, a group of enzymes involved in protein glycosylation. The encoded protein is targeted to Golgi membranes but may be proteolytically processed and secreted. The gene product may also be involved in the increased expression of sialyl Lewis X antigen seen in inflammatory responses. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]