

## Product datasheet for RC232520

### PRKACB (NM\_001242862) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PRKACB (NM_001242862) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PRKACB
Synonyms:	CAFD2; PKA C-beta; PKACB
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC232520 representing NM_001242862 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGCTGCTACTATCTAGTTCTATAAGCTTATATAAAGACAGAAATGAAGCAAGACTGATTCTTCACAGA  
ATAATGCCGGACTGAAGATTTTAAAGGAAAAAACCCTTGGAACAGGTTTCATTTGGAAGATCATGTT  
GGTAAACACAAAGCCACTGAACAGTATTATGCCATGAAGATCTTAGATAAGCAGAAGGTTGTTAACTG  
AAGCAAAATAGAGCATACTTTGAATGAGAAAAGAATATTACAGGCAGTGAATTTTCCTTTCCTTGTTCGAC  
TGGAGTATGCTTTAAGGATAATTCTAATTTATACATGGTTATGGAATATGTCCCTGGGGTGAAATGTT  
TTCACATCTAAGAAGAATTGGAAGGTTTCAGTGAAGCCCATGCACGGTTCTATGCAGCTCAGATAGTGCTA  
ACATTCGAGTACCTCCATTCAGTACCTCATCTACAGAGATCTAAAACCTGAAAATCTCTTAATTGACC  
ATCAAGGCTATATCCAGGTCACAGACTTTGGGTTTGCCAAAAGAGTTAAAGGCAGAACTGGACATTATG  
TGGAACTCCAGAGTATTTGGCTCCAGAAAATAATCTCAGCAAGGGCTACAATAAGGCAGTGGATTGGTGG  
GCATTAGGAGTGCTAATCTATGAAATGGCAGCTGGCTATCCCCATTCTTTCAGACCAACCAATTCAGA  
TTTATGAAAAGATTGTTTCTGAAAGGTCGGATTTCCATCCCCTTCAGTTCAGATCTCAAGGACCTTCT  
ACGGAACCTGCTGCAGGTGGATTTGACCAAGAGATTTGAAAATCTAAAGAATGGTGTGAGTATAAAA  
ACTCACAAGTGGTTTGCCACGACAGATTTGATTGCTATTTACCAGAGGAAGGTTGAAGCTCCATTCATAC  
CAAAGTTTAGAGGCTCTGGAGATACCAGCAACTTTGATGACTATGAAGAAGAAGATATCCGTGTCTCTAT  
AACAGAAAAATGTGCAAAAGAATTTGGTGAATTT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC232520 representing NM\_001242862  
Red=Cloning site Green=Tags(s)


MLLLSSSISLYKDRNEARLISSQNNAGLEDFERKKTGLTGSFGRVMLVKHKATEQYYAMKILDKQKVVKL  
 KQIEHTLNEKRILQAVNFPFLVRLEYAFKDNSNLYMVMEYVPGGEMF SHLRRIGRFSEPHARFYAAQIVL  
 TFEYLHSLDLIYRDLKPENLLIDHQGYIQVDFGFAKRVKGRWTLCGTPEYLAPEIILSKGYNKAVDWW  
 ALGVL IYEMAAGYPPFFADQPIQIYEKIVSGKVRFP SHFSSDLKDLLRNLLQVDLTKRFGNLKNGVSDIK  
 THKWFAATDWIAIYQRKVEAPFIPKFRGSGDTSNFDYEEEDIRVSITEKCAKEFGEF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**ACCN:** NM\_001242862

**ORF Size:** 1014 bp


**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001242862.2</a>
<b>RefSeq Size:</b>	4226 bp
<b>RefSeq ORF:</b>	1017 bp
<b>Locus ID:</b>	5567
<b>UniProt ID:</b>	<a href="#">P22694</a>
<b>Cytogenetics:</b>	1p31.1
<b>Protein Families:</b>	Druggable Genome, Protein Kinase
<b>Protein Pathways:</b>	Apoptosis, Calcium signaling pathway, Chemokine signaling pathway, Dilated cardiomyopathy, Gap junction, GnRH signaling pathway, Hedgehog signaling pathway, Insulin signaling pathway, Long-term potentiation, MAPK signaling pathway, Melanogenesis, Olfactory transduction, Oocyte meiosis, Prion diseases, Progesterone-mediated oocyte maturation, Taste transduction, Vascular smooth muscle contraction, Vibrio cholerae infection, Wnt signaling pathway
<b>MW:</b>	39.7 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a member of the serine/threonine protein kinase family. The encoded protein is a catalytic subunit of cAMP (cyclic AMP)-dependent protein kinase, which mediates signalling through cAMP. cAMP signaling is important to a number of processes, including cell proliferation and differentiation. Multiple alternatively spliced transcript variants encoding distinct isoforms have been observed. [provided by RefSeq, Jul 2014]

Product images:



Circular map for RC232520