

## Product datasheet for **RC232537**

### ST3GAL3 (NM\_001270460) Human Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	ST3GAL3 (NM_001270460) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	ST3GAL3
Synonyms:	DEE15; EIEE15; MRT12; SIAT6; ST3GALII; ST3Gal III; ST3GalIII; ST3N
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RC232537 representing NM_001270460 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGGACTCTTGGTATTTGTGCGCAATCTGCTGCTAGCCCTCGCCTCTTTCTGGTACTGGGATTTTGT  
ATTATTCTGCGTGAAGCTACACTTACTCCAGTGGGAGGAGACTCCAATTCAGTGGTCTTTTCCTTTGA  
CTCCGCTGGACAAACACTAGGCTCAGAGTATGATCGGTTGGCTTCTCCTGAATCTGGACTCTAACTG  
TTCTCCAAGCCAGCACCCATGTTCTGGATGACTCCTTTTCGCAAGTGGGCTAGAATCCGGGAGTTCTGTG  
CGCCTTTTGGGATCAAAGGTCAAGACAATCTGATCAAAGCCATCTTGTCAGTCACCAAAGAGTACCGCCT  
GACCCCTGCCTTGGACAGCCTCCGCTGCCGCCGCTGCATCATCGTGGCAATGGAGGCGTTCTTGCCAAC  
AAGTCTCTGGGGTACGAATTGACGACTATGACATTGTGGTGAGACTGAATTCAGCACCAGTCAAAGGCT  
TTGAGAAGGACGTGGGCAGCAAACGACACTGCGCATCACCTACCCCGAGGGCGCCATGCAGCGGCTGA  
GCAGTACGAGCGCGATTCTCTTTGTCTCGCCGGCTTCAAGTGGCAGGACTTTAAGTGGTTGAAATAC  
ATCGTCTACAAGGAGAGAGTGAAGTGCATCGGATGGCTTCTGAAATCTGTGGCCACTCGAGTGCCCAAGG  
AGCCCCCTGAGATTCGAATCCTCAACCCATATTTTCATCCAGGAGGCCGCTTACCCTCATTGGCCTGCC  
CTTCAACAATGGCCTCATGGGCCGGGGAACATCCCTACCTTGGCAGTGTGGCAGTACCATGGCACTA  
CACGGCTGTGACGAGGTGGCAGTCGACGAGATTTGGCTATGACATGAGCACACCCAACGCACCCCTGCACT  
ACTATGAGACCGTTGCGATGGCAGCCATCAAAGAGTCTGGACGCACAATATCCAGCGAGAGAAAGATT  
TCTGCGGAAGCTGGTAAAGCTCGGTCATCACTGATCTAAGCAGTGGCATC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RC232537 representing NM\_001270460  
 Red=Cloning site Green=Tags(s)

MGLLVFVRNLLLALCLFLVLGFLYSAWKLHLLQWEEDSNSVLSFDSAGQTLGSEYDRLGFLNLNLSKL  
 FSKPAPMFLDDSRKWARIREFVPPFGIKGQDNLKAILSVTKEYRLTPALDSLRCRRCIIVGNNGVLAN  
 KSLGSRIDYDIIVRLNSAPVKGFEDVGSKTTLRITYPEGAMQRPEQYERDSLFLAGFKWQDFKWLKY  
 IVYKERSASDGFWKSVA TRVPKEPPEIRILNPYFIQEAFTLIGLPFNNGLMGRGNIPTLGSVAVTMAL  
 HGCDEVAVAGFGYDMSTPNAPLHYYETVRMAAIKESWTHNIQREKEFLRKLVKARVITDLSSGI


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001270460

**ORF Size:** 1032 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001270460.2</a>
<b>RefSeq Size:</b>	2201 bp
<b>RefSeq ORF:</b>	1035 bp
<b>Locus ID:</b>	6487
<b>UniProt ID:</b>	<a href="#">Q11203</a>
<b>Cytogenetics:</b>	1p34.1
<b>Protein Families:</b>	Secreted Protein, Transmembrane
<b>Protein Pathways:</b>	Glycosphingolipid biosynthesis - lacto and neolacto series, Keratan sulfate biosynthesis, Metabolic pathways
<b>MW:</b>	39.4 kDa
<b>Gene Summary:</b>	The protein encoded by this gene is a type II membrane protein that catalyzes the transfer of sialic acid from CMP-sialic acid to galactose-containing substrates. The encoded protein is normally found in the Golgi apparatus but can be proteolytically processed to a soluble form. This protein is a member of glycosyltransferase family 29. Mutations in this gene have been associated with a form of autosomal recessive nonsyndromic cognitive disability as well as infantile epileptic encephalopathy. Multiple transcript variants encoding several different isoforms have been found for this gene. [provided by RefSeq, Jul 2017]