

Product datasheet for **RC232565**

RAGE (AGER) (NM_001206936) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RAGE (AGER) (NM_001206936) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: AGER
Synonyms: RAGE; SCARJ1; sRAGE
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RC232565 representing NM_001206936
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAGCCGGAACAGCAGTTGGAGCCTGGGTGCTGGTCTCAGTCTGTGGGGGAGTAGTAGGTGCTC
AAAACATCACAGCCCGGATTGGCGAGCCACTGGTCTGAAGTGTAAAGGGGCCCCAAGAAACCACCCCA
GCGGCTGGAATGGAACTGAACACAGGCCGGACAGAAGCTTGAAGGTCTGTCTCCCCAGGGAGGAGGC
CCCTGGGACAGTGTGGCTCGTGTCTTCCCAACGGCTCCCTCTTCCCTCCGGCTGTCTGGGATCCAGGATG
AGGGGATTTCCGGTGCCAGGCAATGAACAGGAATGAAAAGGAGACCAAGTCCAACCTACCGAGTCCGTGT
CTACCAGATTCCTGGGAAGCCAGAAATTGTAGATTCTGCCTCTGAACTCACGGCTGGTGTCCCAATAAG
GTGGGGACATGTGTGTCAGAGGGAAGCTACCCTGCAGGGACTCTTAGCTGGCACTTGGATGGGAAGCCCC
TGGTGCCTAATGAGAAGGGAGTATCTGTGAAGGAACAGACCAGGAGACACCCTGAGACAGGGCTCTTAC
ACTGCAGTCGGAGCTAATGGTGACCCAGCCCGGGGAGGAGATCCCCGTCCCACCTTCTCCTGTAGCTTC
AGCCCAGGCCTTCCCGACACCGGGCCTTGCGCACAGCCCCATCCAGCCCCGTGTCTGGGAGCCTGTGC
CTCTGGAGGAGGTCCAATTGGTGGTGGAGCCAGAAGGTGGAGCAGTAGCTCCTGGTGAACCGTAACCCT
GACCTGTGAAGTCCCTGCCAGCCCTCTCCTCAAATCCACTGGATGAAGGATGGCCTCAGGACCAGGGAA
CCTACAGCTGTGTGGCCACCCATTCCAGCCACGGGCCCCAGGAAAGCCGTGCTGTGAGCATCAGCATCAT
CGAACCAGGCGAGGAGGGGCCAACTGCAGGTGAGGGGTTTGATAAAGTCAGGGAAGCAGAAGATAGCCCC
CAACACATGTGACTGGGGGATGGTCAACAAGAAAGGAATGGAAGGCCCCAGAAAACCAGGAGGAAGAGG
AGGAGCGTGCAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RC232565 representing NM_001206936
 Red=Cloning site Green=Tags(s)

MAAGTAVGAWVLVLSLWGAVVGAQNITARIGEPLVLKCKGAPKKPPQRLEWKLNTGRTEAWKVLSPQGGG
 PWDSVARVLPNGSLFLPAVGIQDEGIFRCQAMNRRNGKTKSNYRVRVYQIPGKPEIVDSASELTAGVPNK
 VGTCVSEGSYPAGTLSWHLDGKPLVPNEKGVSVKEQTRRHPETGLFTLQSELMVTPARGGDRPTFSCSF
 SPGLPRHRALRTAPIQPRVWEPVPLEEVQLVVEPEGGAVPGGTVTLTCEVPAQPSQIHWKMDGLRTRE
 PTAVWPPIPATGPRKAVLSASASSNQARRGQLQVRGLIKSGKQKIAPNTCDWGDGQQRNGRPPQKTRRKR
 RSVQN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001206936

ORF Size:	1065 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001206936.1 , NP_001193865.1
RefSeq Size:	1411 bp
RefSeq ORF:	1068 bp
Locus ID:	177
UniProt ID:	Q15109
Cytogenetics:	6p21.32
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	38.8 kDa
Gene Summary:	The advanced glycosylation end product (AGE) receptor encoded by this gene is a member of the immunoglobulin superfamily of cell surface receptors. It is a multiligand receptor, and besides AGE, interacts with other molecules implicated in homeostasis, development, and inflammation, and certain diseases, such as diabetes and Alzheimer's disease. Many alternatively spliced transcript variants encoding different isoforms, as well as non-protein-coding variants, have been described for this gene (PMID:18089847). [provided by RefSeq, May 2011]